



COMPANY WITH QUALITY SYSTEM CERTIFIED BY DNV = ISO 13485 =

CARMEX

Mobile C arm unit for fluoroscopy and radiography



PRODUCT DATA

Product data CARMEX – Rev. 01 (January 2018)

System Description3Technical Features9Classifications14Installation and warranty14

DESCRIPTION

CARMEX is a mobile X-ray device with image intensifier manufactured for use in orthopaedic, traumatology, abdominal, urology operating theatres.

CARMEX is composed by a Stand and a Monitor Unit, and depending on user needs it can be provided with standard Video Memory (video memory for acquiring and displaying images with last image hold on 2 monitors) or with optional Workstation with Video Processor (video processor: for image aquisition, display and post-processing).

The equipment allows for acquisition of X-ray images in the following modes:

- Continuous fluoroscopy;
- Pulsed fluoroscopy;
- Digital Radiography;
- Digital Subtraction Angiography (D.S.A.), for VIDEO PROCESSOR version only, (*);
- Standard radiography with CR and/or cassette-film.

CARMEX basic version:



CARMEX C-Arm is extremely intuitive to use and easy to move, even in very difficult and restricted areas thanks to its compact and light mechanical structure (extensive use of aluminium alloy). The reduced footprint helps the unit to move around in very small environments.

CARMEX C-Arm is perfectly balanced and does not require continuous use of brakes in normal use. Wheels are equipped with cable pushers and controlled/braked with a single handle.



Thanks to the particular unit balance, the front wheels can be easily lifted in order to overcome any small obstacle that system encounters during its movement.

A touch screen console, placed in the top part of the system, gives the operator the complete control of all operative parameters settings and viewing all messages of alarm/fault conditions.



The X-ray generator is microprocessor controlled and operates at High Frequency (40kHz) to guarantee precise and stable high tension generation. Depending on the monoblock anode type (either stationary or rotating) two versions of CARMEX is available as "F" version with fixed anode and the "R" version with rotating anode. Standard generator power is 3,5 kW (for both "F" and "R" versions) and for "R" version an optional 5 kW is available for specific applications as spine surgery.

Thanks to the Real Time processing tools like automatic Brightness/Contrast, Edge Enhancement and Recursive Filtering for noise reduction and thanks to the 1k x 1k CCD sensor, Image quality is astonishing.

Product Data



Available accessories

STANDARD

- Operativity keyboard for patient data information entry
- Multifunction dual fluoroscopy footswitch
- Radiography hand switch with extension cable
- Set of covers for the entire system with the possibility of sterilization
- Remote control for reviewing images

OPTIONAL

- Laser Localizer
- Dosimeter (DAP Chamber) with presentation on final image in case of digital workstation
- Automatic Contrast/Brightness in area's with metal implants
- Thermal Printer
- 24x30 Cassette holder

OPTIONAL DICOM FEATURES

- DICOM Store
- DICOM Print
- DICOM Worklist
- DICOM Media Interchange (burn data to CD/DVD)
- DICOM Modality Performed Procedure Step (MPPS)
- DICOM Query/Retrieve
- DICOM Storage/Commitment.

TECHNICAL FEATURES

		MECHAN	IICAL DATA (C- ARM)		
		"F" version		"R" version	
I.I. size		9"	12"	9"	12"
Vertical rur	ı (motorized)		45	cm	
Horizontal	run	22 cm			
Wig-wag		<u>+</u> 12°	<u>+</u> 10°	<u>+</u> 12°	±10°
Lateral arc rotation		±270° (total 540°)			
Orbital arc rotation		130° (90°+ 40°)			
S.I.D:		96.9 cm	95.5 cm	96.9 cm	94.2 cm
Arc depth		71,5 cm	68,5 cm	71,5 cm	68,5 cm
Free space between x-ray tube head and i.i.		78,3 cm	73,3 cm	78,3 cm	73,3 cm
Min/max le	ngth	1790/2425 mm			
Max width		775 mm			
Min/max height		1785/2270 mm	1852/2302 mm	1785/2270 mm	1852/2302 mm
	C-arm stand	273 kg	282 kg	280 kg	290 kg
Weight	Monitor unit	115 kg			
Movement		Manual – Parallel steering back wheels with parking position for movements along the patient table			
		Rotating anterior wheel			

ENVIRONMENTAL CONDITIONS (DELIVERY AND STORAGE CONDITIONS)		
Temperature	-10 ÷ 55° C (suggested: 0°÷ 40°C)	
Relative humidity	20 ÷ 90%	
Atmospheric pressure	500 hPa ÷ 1060 hPa	

ENVIRONMENTAL CONDITIONS (OPERATING CONDITIONS)		
Temperature	10 ÷ 40° C	
Relative humidity	30 ÷ 70%	
Atmospheric pressure	700 hPa ÷ 1060 hPa	

	"F" version	"R" version	
GENERATOR			
Power	3,5 kW	3,5 kW (5 kW optional)	
Туре	Microprocess	or controlled	
Generator operating frequency	40	kHz	
kV range	40 ÷ 110 kV	40 ÷ 120 kV	
Max. Radiographic current	6	5 mA	
Max. Image Tube current (fluoro mode)	6 mA (continuous)	
	8 m/	(pulsed)	
Max mA in digital radiography	12 mA	14 mA	
mAs range	1 ÷ 250 mAs		
Anatomical Programs (APR)	Yes, 20 available anatomical programs. More than 1000 combinations.		
RX MONOBLOCK			
Туре	Fixed anode	Rotating anode	
Speed	-	3000 rpm @ 50 Hz	
Focus	0,5x0,5 mm – 1,5x1,5 mm	0,3x0,3 mm – 0,6x0,6 mm	
Anode heat capacity	52 kHU	200 kHU	
Anode continuous dissipation	660 HU/sec	1 kHU/sec	
Monoblock heat capacity	685 kHU (500kJ)	800 kHU (585kJ)	
Additional filtration	1 mm Al		
Total filtration	>3 mm Al		

COLLIMATOR		
Type: Aperture: Control:	From stand control panel	

IMAGE SYSTEM

IMAGE INTENSIFIER	9"	12"	
Model	E 5830SD	TH 9432 HP2	
Image field:			
Nomina	l: 9"/6"/4"	12"/9"/6"	
Workin	g: 215 mm/160 mm/120 mm	290 mm/215 mm/160 mm	
Output image diameter	20 mm	25 mm	
Typical central resolution	64, 60, 52 lp/ cm (4", 6", 9")	56, 50, 44 lp/cm (6", 9", 12")	
Contrast ratio	30:1, 25:1, 19:1 (4", 6", 9")	30:1, 25:1, 22:1 (6", 9", 12")	
DQE @59,5 keV	65% (IEC standard)		
Scintillator	Caesium Io	Caesium Iodide (Csl)	
Conversion factor (typical)	250 Cd/n	250 Cd/m² / mR/s	

GRID	
SID	90 cm
Ratio	8:1
Line rate	103 lines/inch (40 l/cm)

TV MONITOR	N°2 Medical Grade Color 19" LCD	
Technology	TFT	
Angle of view	Vertical and horizontal 178°	
N° of lines	1024	
Resolution	1280x1024	
Contrast ratio 800:1 (typ.)		
Brightness	(monochrome monitor 800 Cd/m ²⁾ 400 Cd/m ²	
Size and dimensions	432x467x251 mm	
Max consumption	80 W	
Height adjustment	20 cm	

TV CAMERA	Standard	High Res (Opt)
Resolution (pixel)	980x980 pixel	1024x1024 pixel
Standard video	Digital serialized	Digital serialized
A/D Conversion	10 bit	12 bit
Bandwidth	20 MHz	40 MHz
Scansion	1:1	1:1
Resolution		2.2 lp/mm

	MEMORY/VIDEO PROCESSOR	
	Video memory	Video processor
Hard disk memory	L.I.H. 30.000 images memory 256 GB capacity Hard-disk	1TByte (average of 360.000 images)
Acquisition frame rate	1 frame/sec (pulsed fluoroscopy)	30 frame/sec
	x-ray time: 880 ms (one-shot fluoroscopy)	
Operative System	-	Windows Embedded
REAL TIME PROCESSING		
Max Opacification	-	X (*)
Virtual collimator	Х	Х
Contrast / Brightness control	-	Х
Gray scale inversion	_	Х
Edge enhancement	Х	Х
Subtraction (D.S.A.)	_	X (*)
Noise reduction	Х	Х
Recursive Filter	Х	Х
Horizontal image flip	X	Х
Digital rotation	Х	Х
Virtual/Digital Collimation	-	Х
Last Image Hold	Х	Х
POST PROCESSING		
Patient & Exam Data Management	Х	Х
Horizontal image flip	Х	Х
Contrast / Brightness control	-	Х
Edge enhancement	-	Х
Reverse polarity (grey scale)	-	Х
Fluoro loop	-	Х
Digital Zoom	-	Х
360° Digital Rotation		
Subtractions (D.S.A.)	_	X (*)
Landmarking	_	X (*)
Digital zoom (2x)	_	Х
Pixel shift	-	X (*)
Text annotations, length, angles	-	Х
Overview (mosaic)	-	Х
Print out	X (*)	X (*)
USB export		Х

MEMORY/VIDEO PROCESSOR		
	Video memory	Video processor
DICOM FUNCTIONS		
Store SCU	-	X (*)
Print SCU	-	X (*)
Worklist SCU	-	X (*)
CD/DVD Recorder w DICOM viewer	-	X (*)
USB Recorder w DICOM viewer	-	X (*)
Modality Performed Procedure Step	-	X (*)
Query/Retrieve	-	X (*)
Storage/Commitment	-	X (*)
DICOM transmission type		Cable or Wireless (*)
Printer	Analogical / digital Sony printer models UP970 and UP 990	Analogical / digital Sony printer models UP970 and UP 990

(*) Optional

CERTIFICATIONS

The CARMEX mobile x-ray unit is a II b class device, according to 93/42 CEE European directive. The device is made according to directives 93/42/EEC, the harmonized norms applicable are:

CEI EN 60601-1	CEI EN 60601-1-6
CEI EN 60601-1-2	CEI EN 60601-2-28
CEI EN 60601-1-3	CEI EN 60601-2-54

INSTALLATION AND WARRANTY

CARMEX mobile x-ray units need only the installation of the monitors up to the monitor trolley, the connection of the two parts (C-arm and monitor unit) and the power supply connection with single phase 230 Vac, earthen.

All ITALRAY device has one-year factory warranty, from delivery date. The customer can extend the warranty by subscription of technical service agreement.

ITALRAY reserves the right to make modifications without any prior notice.



