



Schaerer Configuration and Sales Manual






Axis 500-700 S1
2022



CSF - Carbon Spine Frame



MTS - Modular Traction System

Excellence Through Sustainable Experience



Model X
1912



Nova X-Ray
1957



S-200
1984



Axis 400-800
1994



Arcus 501-701
2009



Axis 350
2019

Arcus - Premium Line

Arcus 501



Arcus 601



Arcus 701



Axis - Top Line

Axis 500



Axis 600



Axis 700



Axis 350 - Advanced Line

Axis 350 Low



Axis 350 High



OP 830 / 850 - Basic Line

EOS 120



EOS 130



EOS 150



Akrus - Basic Chair Table

SC 5010 ES



SC 5010 HS



SC 5010 SEK



AK 5010 MBS



Operating Tables

	ARCUS 501-701	AXIS 500-700	Axis 350	OP830 / OP850
Technical Specifications	6 – 9	40 – 42	70 – 73	94 – 97
Tables and Components	10 – 11	43 – 44	74 – 75	98 – 100
Accessories for General Surgery	12	45	76	101
Ophthalmology, ENT	13	46	77	101
Abdominal	14	47	78	101
Bariatric, Laparoscopic	15	48	79	101
Nephrectomy, Thorax	16	49	80	101
Rectal, Andrews, Genupectoral	17	50	81	101
Gynaecology / Urology	18	51	82	101
Cardiac	19 – 21	52	–	–
Spine	22 – 24	53 – 55	83	102
Neuro	25 – 26	56 – 57	84 – 85	–
Ortho and Trauma – MTS System	27 – 34	58 – 65	–	–
Ortho and Trauma – ETS System	–	–	86 – 89	–
Ortho and Trauma – Hand / Humerus / Knee	35	66	90	–
Ortho and Traumatology – Shoulder	36	67	91	–
Benefits	37 – 38	68	92	–

Surgical Chairs

	Pages
Akrus – SC 5010 ES / SC 5010 HS / SC 5010 SEK	104 – 108
Akrus – AK 5010 MBS	109 – 110
Akrus – Surgeon's Chair AK480	111
Akrus – GenuFix electric knee clamp	112

Accessories

	Pages
Abdominal Rest	132
Accessories Trolley	125
Anaesthesia Screens	118 – 119
Arm and Leg Guards	116
Arm Holders	117
Arm Rests	115 – 116
Carbon Spine Frame (CSF)	135
Carbon-Fibre Plates and accessories	133 – 134
Clamps	114
Collecting Tray TUR	131
Cushions	148 – 150
DORO® Headrest Systems	137 – 138
ETS - Traction System	142
Gel Pads	151 – 157
Gynaecology handle	132
Hand Surgery Tables	143
Head Rest Systems	127 – 128
Headrest Adapters	128 / 137
Hip Prosthesis Positioning Device	145
Humerus / Traction Device	145

	Pages
Infusion Poles	125
Instrument tray	119
Knee Arthroscopy Supports	144
Knee Cradles	132
Knee Positioner	145
Lateral Supports	123 – 124
Leg Holders «Goepel»	130
Leg Holders adjustable	129
MTS - Modular Traction System	139 – 141
Schure Spine Frame	136
Shoulder Plates	146 – 147
Shoulder Positioning Accessories	147
Shoulder Rests	143
Side Rail Extension	125
Straps	120 – 122
Table Extensions	117
Thigh holder	119
Torso Pads	136
X-Ray Cassette Holder	126
X-Ray Cassette Telescopic Slide	126

Arcus 501 – 701 **Premium Line**

Mobile Operating Tables for all weight
categories and surgical applications
with dedicated Accessoires



Arcus – Technical Specifications

Arcus - Technical Specifications	Arcus 501	Arcus 601	Arcus 701
Table dimensions incl. head-and leg plates with side rails (mm)	L 2120 x W 580	L 2120 x W 580	L 2120 x W 580
X-Ray window (mm)	1240	1240	1240
Height adjustment with pad (mm)	595 - 1200	595 - 1200	595 - 1200
Height adjustment without pad (mm)	545 - 1150	545 - 1150	545 - 1150
Head-down Trendelenburg position	30°	30°	30°
Leg-down/Reverse-Trendelenburg position	30°	30°	30°
Seat plate up (motorized)	65°	65°	65°
Seat plate down (motorized)	40°	40°	40°
Leg plate up (manual) (Arcus 601 motorized)	10°	*80°	-
Leg plate down (manual) (Arcus 601 motorized)	90°	*100°	-
Head plate up (manual)	40°	40°	40°
Head plate down (manual)	45°	45°	45°
Lateral tilting to both sides	20°	20°	20°
Side rail dimensions (mm)	25 x 10	25 x 10	25 x 10
Dimension X-ray cassette channel (mm)	W 373.6 x H 32	W 373.6 x H 32	W 373.6 x H 32
Dimension X-ray tray / max. size cassetts (mm)	L 472 x W 365	L 472 x W 365	L 472 x W 365
Table weight without head and leg plates (kg)	300	310	290
Max. patient weight capacity (kg)	500	500	500
Dynamic weight capacity, all positions (kg)	360	360	360
Pad thickness (mm)	50	50	50
Wheels dimensions (mm)	105	105	105
IPX-Protection against entry of liquids	IPX4	IPX4	IPX4




Segment Movements	Arcus 501	Arcus 601	Arcus 701
Height adjustment up / down	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Trendelenburg / Reverse trendelenburg	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Lateral tilt left / right	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Seat plate up / down	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Leg plate up / down	Manual	Electro-hydraulic In A-position	Manual
Head plate up / down	Manual	Electro-hydraulic In B-position	Manual
Flex / Reflex	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Zero position	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Leg plates swivelling out 90°	Manual	Manual	Manual




Electrical Data




- Electrical mains connection 100-120 / 200-240 V AC, 50-60 HZ
- Control voltage 24 V DC
- Main system: (2 x 12 V DC) capacity 12 AH
- Emergency system: 24 V DC (2 x 12 V DC) capacity 7 AH
- 4 maintenance-free batteries
- Battery charging time approx. 6 to 8 hours
- Battery capacity approx. 20 operations
- Charge indication on the hand control
- Low battery acoustic signal






Arcus – Trendelenburg / Reverse Trendelenburg

Degrees	Highest Position	Lowest Position
 + 30° - 30°	 1905 mm 1575 mm 1325 mm 1085 mm 935 mm	 1300 mm 970 mm 720 mm 480 mm 330 mm




Arcus – Head Rest	Arcus – Leg Plates	Arcus – Beach-Chair-Position
 + 45° - 45°	 + 5° - 90°	 + 65° - 40°

Arcus – Flex	Arcus – Reflex	Arcus – Lateral Tilt
 220°	 115°	 20° 20°

Arcus 601 – With Second Motorized Joint


Arcus 601 – Flex	Arcus 601 – Head or Leg Plates	Arcus 601 – Head or Leg Plates
 + 280°	 - 100°	 100° 80°

Dimensions




Foot dimensions


Arcus 501 – Segments








Arcus 601 – Segments



A - Position



B - Position




Arcus - Range



Without pads:	545 - 1150 mm
With pads:	595 - 1200 mm
Column range	605 mm
Height pads:	50 mm

Arcus – Patient Weight Capacity



Static 500 kg:


- ✓ Standard patient position

Dynamic 360 kg:


- ✓ Patient weight supported in **all positions** and execution of **all movements** of the tabletop

Arcus – X-Ray Window


Standard Configuration




Adding Intermediate Plate (Item: 92231)




With short Carbon Plate (Item: 90612)



With large Carbon Plate (Item: 90597)



Arcus - Controls




This command with cable allows to realize all the movements of the table. In addition it has a function that indicates all the errors that can occur by a system of codes that allows a fast and effective detection.

The emergency control allows all movements to be carried out with a second system completely separated from the central system. This system activates the hydraulic system in case the control stops working or when the batteries are exhausted.

Foot switch (optional)

Item: **92033**

Foot switch for **Arcus**



- 1 - Button column up
- 2 - Button column down
- 3 - Button trendelenburg
- 4 - Button rev. trendelenburg
- 5 - Button tilt right
- 6 - Button tilt left

Electrical drive unit (optional)



Fast Speed

Pressing two buttons at one time initiates the fast speed

Slow Speed

Pressing one direction button initiates the slow speed for manoeuvres in narrow areas and keeping one hand free, or carrying additional items (i.e. infusion stands on wheels).

Arcus - Models

Arcus 501 CrNi

Item: **92226**

Arcus 501 CrNi

Item: **92226D**

Arcus 501 CrNi with drive unit

Item: **91260**

Packaging



- Removable seat plate
- Can be used for all surgeries
- All accessories CSF / MTS or any compatible Schaerer traction device can be attached
- One electric-hydraulically driven tabletop movement for seat plate

Arcus 601 CrNi

Item: **92227**

Arcus 601 CrNi

Item: **92227D**

Arcus 601 CrNi with drive unit

Item: **91260**

Packaging



- Seat plate not removable
- Ideally for Gyn/Obs, Urology and Neurosurgery
- Not for Ortho/Trauma
- Two electro-hydraulically driven table top movements for seat and leg plates

Arcus 701 CrNi

Item: **92228**

Arcus 701 CrNi

Item: **92228D**

Arcus 701 CrNi with drive unit

Item: **91260**

Packaging



- Table same like arcus 501 but without seat plate
- Can be used for all surgeries
- All accessories CSF / MTS or any compatible Schaerer traction device can be attached
- One electric-hydraulically driven table top movement for seat plate.

Arcus - Positions



Maximum free space for C-arm utilization in the basic position



Reverse position with carbon plate for 1550 mm X-ray window (1350 mm 3D)



Gynaecological, urological and laparoscopic interventions
patient weight capacity 360 kg














Lateral positioning for nephrectomy, XLIF, dynamic
patient weight capacity 360 kg




Aluminium or radiolucent headrest system with 3-pin rigid cranial fixation. For supine, prone, lateral and sitting positions



Spine intervention with four post positioning system

Head Plates		Bariatric Leg Plates
Item: 92210 Head rest standard  300 mm 520 mm Without rails	Item: 92211 Head rest double joint  300 mm 520 mm Without rails	Item: 92450 Laterally sliding Leg-plates Arcus 
Leg Plates		
Item: 92053 Leg plates 2-part  700 mm 520 mm Without rails	Item: 92049 Leg plate one-piece rigid  670 mm 520 mm Without rails	Item: 92050 Leg plate one-piece  700 mm 520 mm Without rails
Intermediate Plates		
Item: 92231 Intermediate plate long  450 mm 520 mm Without rails	Item: 92232 Intermediate plate short  230 mm 520 mm Without rails	Item: 92051 Seat bar short  280 mm 520 mm Without rails
Carbon Plates		
Item: 90597 Carbon-fibre plate 1400mm for Arcus 501/701 (Instead of Seat plate)  1400 mm 130 mm 520 mm	Item: 90596 Carbon-fibre plate 850 mm for Arcus 501/701 (Instead of Seat plate)  850 mm 130 mm 520 mm	Item: 90612 Carbon-fibre plate 850 mm for Arcus 501/601 (Instead of Leg plates)  850 mm 130 mm 520 mm

Basic Configuration





No	Item	Qty	Description
1	92226 92227	1	Arcus 501 CrNi (or) Arcus 601 CrNi
2	92210 92211	1	Head rest standard (or) Head rest double joint
3	92053 92050 92049	1	Leg plates 2-part (or) Leg plate one-piece (or) Leg plate one-piece, rigid
	91260 91261	1 1	Packaging OT-Table WITHOUT ramp Reusable exit ramp
Optional Accessories			
4	92231	1	Intermediate plate, long
5	92232	1	Intermediate plate, short

Accessories for General Surgery



No	Item	Qty	Description
6	90321	1	Anaesthesia screen, rigid (+)
	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
6	90323	1	Anaesthesia screen with central articulation (+)
	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
6	90635	1	Anaesthesia screen, with extendable rods (+)
	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
6	90636	1	Anaesthesia screen central art., ext. rods (+)
	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
7	90520	2	Arm rest (or)
			4646.05 Arm support quick central fixation (+)
			5.890-VE Fastening straps (1 pair)
8	800-0057	1	Restraint Strap
9	92021	1	Arm strap with fixing screw
10	90500	1	Infusion pole with clamp, height-adjustable
11	90332	1	Instrument tray with clamp, height-adjustable
12	90360	1	Accessories trolley (+)
	91300	1	Packaging for Accessories trolley
Dedicated Accessories			
13	92037	1	Universal Adapter for DORO® basic unit
14	92027	1	Collecting tray TUR for Arcus 501/601
15	90260	1	Hand surgery table with fixing screw
16	92006	1	X-ray cassette frame for arcus

Ophthalmology / ENT – Configuration A



Ophthalmology / ENT – Configuration B



No	Item	Qty	Description
1	92027	1	Adapter for DORO® basic unit, for Arcus
2	90045	1	Horseshoe-headrest holder 3 joints
3	1002.72A0	1	Horseshoe-shaped head rest, two pieces
4	1004.75A0	1	Wrist hand support for Ophthalmology

Abdominal – Configuration A

1



Item: 90520

Item: 4646.05
Item: 5.890-VE

2



Item: 90325

Item: 90506

Item: 90310

3



Item: 90330

Item: 90506

Item: 90310



No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90325 90506	1 1	Arm strap with fixing screw Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	90330 90506	2 2	Leg holder "Goepel" (1 piece) (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp

Abdominal – Configuration B

1



Item: 90520

Item: 4646.05
Item: 5.890-VE

2



Item: 90325


Item: 90506

Item: 90310

3





Adjustable Leg holder



No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90325 90506	1 1	Arm strap with fixing screw Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	800-0342-PR 800-0338-EU	1 2	Adjustable Leg holder for patients up to 182kg (or) 800-0342-PL for patients up to 272 kg (or) 800-0364 for patients up to 454 kg (+) Springloc Clamp (order 2 pieces)


Bariatric, Laparoscopic – Configuration A

1  Item: 90520 (or)  Item: 4646.05 + 5.890-VE	2  Item: 90079 (+)  Item: 81501	3  Item: 800-0055 (+)  Item: 800-0338-EU	4  Item: 92038
--	--	--	---



No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90079 81501	1 2	Knee cradles (1 pair), with pads (+) Clamp, Knee cradles (order 2 pieces)
Optional			
3	800-0055 800-0338-EU	1 2	Nissen straps, 272 kg patient w. cap. (1 pair) (+) Springloc Clamp
4	92038	4	Table extension for Arcus (1 piece),

Bariatric, Laparoscopic – Configuration B

1  Item: 90520  Item: 4646.05 + 5.890-VE	2  (+) Adjustable Leg holder	3  Item: 800-0055 (+)  Item: 800-0338-EU	4  Item: 92038
---	--	--	---



No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	800-0342-PR 800-0338-EU	1 2	Adjustable Leg holder for patients up to 182 kg (or) 800-0342-PL for patients up to 272 kg (or) 800-0364 for patients up to 454 kg (+) Springloc Clamp (order 2 pieces)
3	800-0055 800-0338-EU	1 2	Nissen straps, 272 kg patient w. cap. (1 pair) (+) Springloc Clamp
4	92038	4	Table extension for Arcus (1 piece),

Nephrectomy, Thorax

1



Item: **90520**

(or)



Item: **4646.05**
Item: **5.890-VE**



2



Item: **90325**

(+)



Item: **90506**
(or)



Item: **90310**

3



Item: **90351**

(+)



Item: **90350**

(+) **(or)**



Item: **90539**



Item: **90314**

4



Item: **90352**

(+)



Item: **90350**

(+) **(or)**



Item: **90539**



Item: **90314**

5



Item: **90609**

(+)



Item: **90350**



(+) **(or)**



Item: **90539**




Item: **90314**



No	Item	Qty	Description
1	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90325 90506	1 1	Arm rest high, without clamp (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	90351 90350 90539	1 1 1	Side and back rest, with pad (+) Universal holder for side and back rest (+) Radial clamp with fixing screw, rotating (or) 90314 Rapid-action radial setting clamp
4	90352 90350 90539	1 1 1	Side and back supporting roll Universal holder for side and back rest (+) Radial clamp with fixing screw, rotating (or) 90314 Rapid-action radial setting clamp
5	90609 90350 90539	1 1 1	Sacrum/Pubis/Sternum support, with pad (+) Universal holder for side and back rest (+) Radial clamp with fixing screw, rotating (or) 90314 Rapid-action radial setting clamp

Rectal, Andrews, Genupectoral

1  (or)  Item: 90520 Item: 4646.05 + 5.890-VE	2  (+)  Item: 90506 (or)  Item: 90310
3  (+)  Item: 90079 Item: 81501	4  (or)  Item: AP022 Item: QPHPD-ADULT 1



No	Item	Qty	Description
1	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90277 90506	1 2	Abdominal rest, with pad, adjustable (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	90079 81501	1 2	Knee cradles (1 pair), with pads (+) Clamp, Knee cradles (order 2 pieces)
4	AP022	1	Gel Headrest for prone position (or) QPHPD-ADULT 1 – Head rest (single use)

Gynaecology / Urology – Configuration A

1



Item: 90330



Item: 90506

(or)

Item: 90310

2



Item: 92037

3



Item: 90520

(or)



Item: 4646.05 + 5.890-VE



No	Item	Qty	Description
1	90330	2	Leg holder "Goepel" (1 piece) (+)
	90506	2	Radial setting clamp with fixing screw (or)
	90310		Rapid-action radial setting clamp
2	92037	1	Collecting tray TUR for Arcus (or)
3	90520	2	Arm rest (or)
	4646.05		Arm support quick central fixation (+)
	5.890-VE		Fastening straps (1 pair)

Gynaecology / Urology – Configuration B

1



Adjustable Leg holder

2



Item: 92037

3



Item: 90520

(or)



Item: 4646.05 + 5.890-VE

4



Item: 800-0074-S



No	Item	Qty	Description
1	800-0342-PR	1	Adjustable Leg holder for patients up to 182 kg (or)
	800-0342-PL		for patients up to 272 kg (or)
	800-0364		for patients up to 454 kg (+)
	800-0338-EU	2	Springloc Clamp (order 2 pieces)
2	92037	1	Collecting tray TUR for Arcus
3	90520	2	Arm rest with clamp (or)
	4646.05		Arm support quick central fixation (+)
	5.890-VE		Fastening straps (1 pair)
Optional			
4	800-0074-S	1	Great White® Stirrup Dolly

CFP Carbon Fibre Plates

Item: 90597


Carbon-fibre plate 1400 mm
for Arcus 501/701

Patient weight capacity: 130kg

Dimensions: 1400 x 520 x 130 mm


(Instead of Seat plate)

(Not for Arcus 601)



Item: 90593

Side rail 25 x 10 mm (1 piece)



Item: 90536


Arm rest carbon-fibre for Arcus






No	Item	Qty	Description
1	92228 92226	1	Arcus 701 CrNi <i>(or)</i> Arcus 501 CrNi <i>(Without seat plate)</i>
2	92210 92211	1	Head rest standard <i>(or)</i> Head rest double joint
3	90597	1	Carbon-fibre plate 1400 mm
4	90593	2	Side rail 25 x 10 mm (1 piece)
5	90536	1	Arm rest carbon-fibre for Arcus with pad
Optional			
6	92053 92050 92049	1	Leg plates 2-part <i>(or)</i> Leg plate one-piece <i>(or)</i> Leg plate one-piece, rigid <i>(or)</i>
	92231	1	Intermediate plate, long <i>(or)</i>
	92232	1	Intermediate plate, short



A - Position



B - Position




CFP Carbon Fibre Plates

Item: 90596

Carbon-fibre plate 850 mm
for Arcus 501/701.

Patient weight capacity: 220kg
Dimensions: 850 x 520 x 130 mm
(Instead of Seat plate)
(Not for Arcus 601)



Item: 90593

Side rail 25 x 10 mm (1 piece)

4



Item: 90536

Arm rest carbon-fibre for Arcus


5






No	Item	Qty	Description
1	92228 92226	1	Arcus 701 CrNi (or) Arcus 501 CrNi (Without seat plate)
2	92210 92211	1	Head rest standard (or) Head rest double joint
3	90596	1	Carbon-fibre plate 850 mm
4	90593	2	Side rail 25 x 10 mm (1 piece)
5	90536	1	Arm rest carbon-fibre for Arcus with pad
Optional			
6	92053 92050 92049	1	Leg plates 2-part (or) Leg plate one-piece (or) Leg plate one-piece, rigid (or)
	92231	1	Intermediate plate, long (or)
	92232	1	Intermediate plate, short



A - Position



B - Position



CFP Carbon Fibre Plates

Item: 90612


Carbon-fibre plate 850 mm
for Arcus 501/601.

Patient weight capacity: 90kg

Dimensions: 850 x 520 x 130 mm

(Instead of the Leg plate)

(Not for Arcus 701)



Item: 90593

Side rail 25 x 10 mm (1 piece)

4



Item: 90536

Arm rest carbon-fibre for Arcus

5






No	Item	Qty	Description
1	92228 92226	1	Arcus 701 CrNi (or) Arcus 501 CrNi (Without seat plate)
2	92210 92211	1	Head rest standard (or) Head rest double joint
3	90612	1	Carbon-fibre plate 850 mm
4	90593	2	Side rail 25 x 10 mm (1 piece)
5	90536	1	Arm rest carbon-fibre for Arcus with pad
Optional			
6	92053 92050 92049	1	Leg plates 2-part (or) Leg plate one-piece (or) Leg plate one-piece, rigid (or)
	92231	1	Intermediate plate, long (or)
	92232	1	Intermediate plate, short



A - Position




B - Position



Spine – Configuration

1  Item: NO 165	2  Item: NA 401	3  Item: 90520 Item: 4646.05 Item: 5.890-VE
4  Item: 92027 Item: 90456 Item: 90041 Item: 90003 Item: QPHPD-ADULT 1 Item: AP022		
5  Item: 90612	6  Item: 800-0235	7  Item: 800-0074-SF



No	Item	Qty	Description
1	NO 165	1	Torso Pad "Onsernone"
2	NA 401	1	Leg support wave cushion 65x50x13cm
3	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
4	92027 90456 90041 90003	1 1 1 1	Adapter for DORO® basic unit, for Arcus (+) Hexagonal adapter with square adapter (+) Head rest holder, with lateral adjustment (+) Head rest (horseshoe) with pad (or) QPHPD-ADULT 1 - QView™ Head positioning device (single use) (or) AP022 Gel Headrest for prone position
Optional Accessories			
5	90612	1	Carbon-fibre plate 850 mm
6	800-0235	1	SchureSpine Frame, Radiolucent
7	800-0074-SF	1	Schure Spine Frame Dolly

CSF Carbon Spine Frame

1



Item: 91600

2



Item: 90520

(or)



Item: 4646.05
Item: 5.890-VE

3



Item: 91603

(or)



Item: 91604

4



Item: AP022

(or)



Item: QPHPD-ADULT 1

Item: 91600


Schaerer Carbon Spine Frame

Patient weight capacity:

Arcus: 240 kg




No floor
support




360°
Radiolucent

No	Item	Qty	Description
1	91600	1	Schaerer Carbon Spine Frame (CSF)
2	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
Optional Accessories			
3	91603 91604	1	Side rail adapter short (or) Side rail adapter long
4	AP022	1	Gel Headrest for prone position (or) QPHPD-ADULT 1 - QViewTM Head positioning device (single use)

CSF - Carbon Spine Frame DORO® Radiolucent Head Rest System






(or)



Radiolucent			
No	Item	Qty	Description
6	1101.021	1	DORO® LUCENT Basic support
7	1101.031	1	DORO® LUCENT Clamp intermediate piece
8	1101.026	1	DORO® LUCENT Intermediate piece
9	3033-00	1	DORO® Swivel adapter, radiolucent
10	3034-00	1	DORO® Skull clamp, radiolucent Quick-Rail



Aluminium			
No	Item	Qty	Description
11	3001-00	1	DORO® Base unit aluminum
12	3002-00	1	DORO® Swivel adaptor aluminum (or)
13	1001.001	1	DORO® Skull clamp QR3 aluminum

Technical Data




Patient weight capacity	240 kg
Overall length	1205 mm
Length 360° radiolucent	805 mm
Frame length	1052 mm
Width frame	480 mm

Neuro – Configuration A - Aluminium Skull Clamp System



No	Item	Qty	Description
1	92027	1	Adapter for DORO® Basic unit, Arcus (or) 3011-00 Rail adapter (for all tables)
2	3001-00	1	DORO® Base unit aluminum
3	3002-00	1	DORO® Swivel adaptor aluminum (or) DORO® Ball adapter, allows 360° rotation
4	1001.001	1	DORO® Skull clamp QR3 aluminum
5	3009-00	1	DORO® Horse shoe head rest for adults
6	3008-00	1	DORO® Horse shoe head rest for children
		1	Pins reusable (or) Pins disposable
Additional accessories for sitting position			
7	3007-00 90506	1	DORO® Rail adapter for sitting position (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp



Fluoroscopy

Pins reusable



Item: 3005-00 (Stainless Steel)
Item: 3005-50 (Titanium)
Reusable skull pins for **adults**
(3 pieces)



Item: 3004-00 (Stainless Steel)
Item: 3004-50 (Titanium)
Reusable skull pins for **children**
(3 pieces)

Pins disposable



Item: 1106.003
DORO® Disposable skull pins, for **adults**,
(36 pieces)



Item: 3006-10
DORO® Disposable skull pins, for **children**,
(36 pieces)

Pins disposable



Item: 3006-20
DORO® Disposable skull pins, for **adults**,
(36 pieces)





Item: 3006-30
DORO® Disposable skull pins, for **children**,
(36 pieces)

Neuro – Configuration B - Radiolucent Skull Clamp System


Item: **1101.020** (Complete Set)

DORO® Radiolucent Head Rest System Compatible with all tables
(No adapter needed)



No	Item	Qty	Description
1	1101.021	1	DORO® LUCENT Basic support
2	1101.031	1	DORO® LUCENT Clamp intermediate piece
3	1101.026	1	DORO® LUCENT Intermediate piece
4	3033-00	1	DORO® Swivel adapter, radiolucent
5	3034-00	1	DORO® Skull clamp, radiolucent with Quick-Rail™
6		1	Pins reusable (or) Pins disposable
Additional accessories for sitting position			
7	3007-00 90506	1	DORO® Rail adapter for sitting position (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
8	3001-00	1	DORO® Base unit aluminum

Sitting Position



Pins reusable



Item: **3005-00** (Stainless Steel)
Item: **3005-50** (Titanium)
Reusable skull pins for **adults**
(3 pieces)



Item: **3004-00** (Stainless Steel)
Item: **3004-50** (Titanium)
Reusable skull pins for **children**
(3 pieces)

Pins disposable



Item: **1106.003**
DORO® Disposable
skull pins, for **adults**,
(36 pieces)



Item: **3006-10**
DORO® Disposable
skull pins, for **children**,
(36 pieces)

Pins disposable



Item: **3006-20**
DORO® Disposable
skull pins, for **adults**,
(36 pieces)



Item: **3006-30**
DORO® Disposable
skull pins, for **children**,
(36 pieces)


MTS Modular Traction System (One Surgical Table - Multiple Disciplines)

Applications


<u>General Surgery</u>	<u>Femur</u>	<u>Tibia</u>	<u>Hip</u>	<u>Lateral Position</u>	<u>Shoulder</u>
					

Dimensions in retracted position (in mm)

- MTS




- MTS + MIS




Dimensions in extended position (in mm)


- MTS




- MTS + MIS









Dimensions seat plate (in mm)




Dimensions MTS-Trolley (in mm)






MTS - Femur			
No	Images	Item	Description
1		90658	Seat Plate ADVANCED, Arcus 501/701
2		90661	Supplement bar with plate, left, Arcus 501/701
3		90662	Supplement bar with plate, right, Arcus 501/701
4		90971	Perineal rod
5		90990	Perineal rod for arthroscopy
6		90958	Offset holder
7		90983	Additional joint
8		90981	Traction bar Metal
9		90974	Traction bar ADVANCED

MTS - Femur			
No	Images	Item	Description
10		90978	Multiclamp (For Traction bar ADVANCED only)
11		90928	Multiclamp slide (For Traction bar ADVANCED only)
12		90980	Traction unit
13		90895	Traction unit, vertical shift (allows vertical movement during traction)
14		90763	Extension shoe MTS, white (1 piece)
15		RO.0162 .2018	Traction shoe with side rail interface
16		RO.0148 .2018	Pad set incl. fixation straps (10 pieces)
17		90330 (+) 90506	Leg holder "Goepel" (1 piece) (+) 90506 Radial setting clamp with fixing screw
18		90994	Extension bow holder, for Kirschner- wire or pin

MTS – Tibia additional Items

No	Images	Item	Description
19		90780	Central rail adapter
20		90897	Radial setting clamp with star screw
21		90759	Holder for Knee rest
22		90896	Adapter for lower leg surgery
23		92699 92698	Foot plate Adult Foot plate for Child
24		90193	Telescopic support



MTS - Shoulder additional Items

No	Images	Item	Description
25		90688	Head rest system, Arcus 501/701
26		F-90023	Head Pad for MTS (10 pieces)
27		90745	Extension plate (for tall patients) Arcus 501/701
28		90677	Lateral ext. plate for head rest system, left, Arcus 501/701
29		90675	Lateral ext. plate for head rest system, right, Arcus 501/701

MTS - Lateral Position additional Items

No	Images	Item	Description
30		92611	ETS Counter Traction Unit for lateral position

Optional Items

31		90670	General Surgery pad set for MTS Seat plate
32		90673	Gyn. cut-out plate, Arcus 501/701

MTS - MIS Items

No	Images	Item	Description
34		90738 91298	Bar MIS for MTS system (1 piece) Packaging MIS
35		90775 91299	MIS Trolley Packaging Trolley

MTS - Trolley

No	Images	Item	Description
33		90957 91304	MTS Trolley Packaging Trolley Measurements: Width: 625 mm Length: 890 mm Height: 1610 mm

MTS – Table Configuration

A



Item: **92228**

B




Item: **92210**

C




Item: **92053**

1




Item: **90658**

2





Item: **90661**

3






Item: **90662**




No	Item	Qty	Description
MTS Table			
A	92228	1	Arcus 701 CrNi (or) Arcus 501 CrNi
B	92210	1	Head rest standard (or) 92211 Head rest double joint
C	92053	1	Leg plates 2-part
Accessories			
1	90658	1	Seat Plate ADVANCED, for Arcus
2	90661	1	Supplement bar with plate, left, for Arcus
3	90662	1	Supplement bar with plate right, for Arcus

Femur - Traction with Metal bars (Additional to MTS Table)


4  Item: 90971	6  Item: 90958	7  Item: 90983	8  Item: 90981	
12  Item: 90980	(or) 13  Item: 90895	14  Item: 90763	(or) 15  Item: RO.0162.2018	(+) 16  Item: RO.0148.2018

1 Leg Configuration



No	Item	Qty	Description
4	90971	1	Perineal rod
6	90958	1	Offset holder
7	90983	1	Additional joint
8	90981	1	Traction bar Metal
12	90980	1	Traction unit (or)
13	90895	1	Traction unit, vertical shift
14	90763	1	Extension shoe MTS, white (1 piece) (or)
15	RO.0162.2018	1	Traction shoe with side rail interface (+)
16	RO.0148.2018	1	Pad set straps (disposable) (10 pieces)
17	90330 90506	1	Leg holder "Goepel" (1 piece) (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp


2 Leg Configuration



No	Item	Qty	Description
4	90971	1	Perineal rod
6	90958	1	Offset holder
7	90983	1	Additional joint
8	90981	2	Traction bar Metal
12	90980	2	Traction unit (or)
13	90895	2	Traction unit, vertical shift
14	90763	2	Extension shoe MTS, white (1 piece) (or)
15	RO.0162.2018	2	Traction shoe with side rail interface (+)
16	RO.0148.2018	1	Pad set straps (disposable) (10 pieces)


Femur - Traction with Carbon bars (Additional to MTS Table)

1 Leg Configuration




No	Item	Qty	Description
4	90971	1	Perineal rod
6	90958	1	Offset holder
7	90983	1	Additional joint
9	90974	1	Traction bar ADVANCED
10	90978	1	Multiclamp (or)
11	90928		Multiclamp slide
12	90980	1	Traction unit (or)
13	90895		Traction unit, vertical shift
14	90763	1	Extension shoe MTS, white (1 piece) (or)
15	RO.0162.2018		Traction shoe with side rail interface (+)
16	RO.0148.2018		Pad set straps (dispos.) (10 pieces)
17	90330 90506	1	Leg holder "Goepel" (1 piece) (+) Radial setting clamp with screw (or) 90310 Rapid-action radial setting clamp


2 Leg Configuration



No	Item	Qty	Description
4	90971	1	Perineal rod
6	90958	1	Offset holder
7	90983	1	Additional joint
9	90974	2	Traction bar ADVANCED
10	90978	2	Multiclamp (or)
11	90928		Multiclamp slide
12	90980	2	Traction unit (or)
13	90895		Traction unit, vertical shift
14	90763	2	Extension shoe MTS, white (1 piece) (or)
15	RO.0162.2018		Traction shoe with side rail interface (+)
16	RO.0148.2018		Pad set straps (dispos.) (10 pieces)

Tibia - in addition to Femur Traction with **Carbon** or **Metal** bars

19	20	21	22	23	24
					
Item: 90780	Item: 90897	Item: 90759	Item: 90896	Item: 90699	Item: 90193



No	Item	Qty	Description
19	90780	1	Central rail adapter
20	90897	1	Radial setting clamp with star screw
21	90759	1	Holder for knee rest
22	90896	1	Adapter for lower leg surgery
23	92699	1	Foot plate / 92698 foot plate for child
24	90193	1	Telescopic support

Lateral Position - in addition to Femur Traction with **Carbon** or **Metal** bars




Item: 92611



Item: 90325 (+) 90506





Item: 90520



No	Item	Qty	Description
30	92611	1	Counter Traction Unit for lateral position
	90325	1	Arm rest high, without clamp (+) 90506
	90520	1	Arm rest (or) 4646.05 (+) 5.890-VE

Shoulder (Additional to MTS Table)

25  Item: 90688	26  Item: F-90023	27  Item: 90745	28  Item: 90677	29  Item: 90675
---	---	---	--	---








No	Item	Qty	Description
25	90688	1	Head rest system, Arcus 501/701
26	F-90023	1	Head Pad for MTS (10 pieces)
27	90745	1	Extension plate, Arcus 501/701
28	90677	1	Lateral ext. plate for head rest system, left, Arcus 501/701
29	90675	1	Lateral ext. plate for head rest system, right, Arcus 501/701

Possible to combine **Metal** with **Carbon**



MTS - MIS

1  Item: 90658	2  Item: 90661	3  Item: 90662	4  Item: 90971	6  Item: 90958
7  Item: 90983	8  Item: 90981	12  Item: 90980	14  Item: 90763	34  Item: 90738



MIS Trolley – Optional (Without accessories)



Item: 90775

MIS Trolley

(+)

Item: 91299

Packaging MIS Trolley

No	Item	Qty	Description
1	90658	1	Seat Plate ADVANCED, for Arcus
2	90661	1	Supplement bar with plate, left, for Arcus
3	90662	1	Supplement bar with plate right, for Arcus
4	90971	1	Perineal rod
6	90958	1	Offset holder
7	90983	1	Additional joint
8	90981	1	Traction bar Metal
12	90980	1	Traction unit (or)
13	90895		Traction unit, vertical shift
14	90763	2	Extension shoe MTS, white (1 piece) (or)
15	RO.0162.2018		Traction shoe with side rail interface (+)
16	RO.0148.2018		Pad set straps (disposable) (10 pieces)
34	90738 91298	1	Bar MIS, for MTS Packaging MIS

Hand Surgery



1



1

(or)



2

(+)



Item: **90506**
(order 2 pc.)

No	Item	Qty	Description
1	90260	1	Hand surgery table, with fixing screw (or)
2	800-0040	1	Hand surgery table, radiolucent
	90506	2	Radial setting clamp

Humerus / Traction Device



1





No	Item	Qty	Description
1	9907038	1	Humerus and Forearm Positioning Device

Knee Arthroscopy



1



(or)



(+)



(or)

Item: **90357** Item: **90356**

Item: **90539** Item: **90314**



2



(+)



(or)

Item: **90330** Item: **90506** Item: **90310**




1



2

No	Item	Qty	Description
1	90357	1	Leg holder arthroscopy (or)
	90356	1	Knee rest (+)
	90539	1	Radial clamp with fixing screw, rotating (or)
	90314	1	Rapid-action radial clamp, rotating
2	90330	1	Leg holder "Goepel" (1 piece) (+)
	90506	1	Radial setting clamp with fixing screw (or)
	90310	1	Rapid-action radial setting clamp

SSP - Schaerer Shoulder Plate - Configuration A


1  Item: 92655	2  Item: 90456	3  Item: 90041	4  tem: 90566 (or) Item: 90003
--	--	---	---



No	Item	Qty	Description
1	92655	1	Shoulder plate (SSP) schaerer® arcus
2	90456	1	Adapter
3	90041	1	Head rest holder, with lateral adjustment
4	90566 90003	1	Helmet, radiolucent (or) Head rest (horseshoe) with pad

SSP - Schaerer Shoulder Plate - Configuration B

1  Item: 92655	2  Item: 92666	3  Item: 90899	4  Item: F-90023
--	--	---	--



No	Item	Qty	Description
1	92655	1	Shoulder plate (SSP) schaerer® arcus
2	92666	1	Head rest holder
3	90899	1	Head rest
4	F-90023	1	Head pad (10 Pieces)


1 - Maximum capacity

High weight capacity (500kg/360kg) means more safety for all kind of patient positions, all kind of operations, all surgical disciplines and a much more stable table. There is no sliding who limits the weight capacity.



2 - Maximum stability

Self-compensating dynamic floor lock system. The table stands on 4 big banana size stamps and offers best stability compared to others in the market. The dynamic hydraulic floor lock system corrects uneven floors.



3 - Maximum height adjustment

Only a large column range from very low to very high with an advanced quality telescopic column design is able to move safe and stable from very low to very high with a hub of 605mm and a column range of 595-1200mm (545-1150 without cushions).




4 - Maximum accessibility

Maximum free area because of eccentric column design. Even in standard configuration the arcus table has a huge x-ray window of 1240mm. Tables without sliding are more stable and safer than tables with sliding. All mechanical and movable parts of a table makes it less safe, can break and reduce weight capacity. During an operation there is no need for sliding. Eccentric column design has technical advantages and less limitations.



5 - Maximum safety

Arcus 501 – 701 tables are the only tables available with a redundant backup system. Consisting out of a complete 2nd system including a separate power supply, batteries, control devices and motor/pump units. A table without an independent 2nd back-up system may have a safety issue if only one system is not working properly.



Emergency Panel on table



Board with dual components in case of failures



2nd battery (in case of emergency)

Main battery

6 - Maximum patient comfort

Accurate and smooth (slow mode) motion for safe movements during surgery. High quality mattress also for long operations and heavy weight patients. A high-quality pad is important for patients with more weight and long operations. If a pad is too thick or soft, it can cause more pressure inquiries to the patient. You cannot move or reposition the patient on the table that easy. Also, for normal seize/weight patients the quality of the pad is of importance and not primarily the thickness.



7. Maximum flexibility

Mobile operating table for all surgical disciplines. Arcus is a universal table, which can be used for all surgical disciplines. It goes very low and very high and is made for high patient weight capacity. With all the unique features surgeons can work and do all treatments with no limitation. Some disciplines needs tables which goes very low (Spine, sitting positions) and other which goes very high (Ortho). The table offers various upgrade options such as motor drive or digital integration.



8 - Maximum quality

Swiss precision wherever you operate since 1892
Manufactured in Switzerland with design and engineering on highest standards. More than 100 years of experience. Due to high Quality, manufacturing product failure rate is low. This means low running costs and low total costs of ownership and long product lifetime. Our products are manufactured in accordance with certified quality standards. ISO 13485:2003.



9 – Slow mode function

Accurate and smooth (slow mode) motion for safe movements during surgery. In critical situations you can have slow and safe movements of the table e.g. in Laparoscopy procedures.



10 – Emergency Panel on the base

The table is covered and draped with towels. Sometimes there is no good access to the column of the table, and it is easier access to the base of the table (head side). Additionally, the surgeons stand on the side of the column and there is no free access for the staff to the emergency panel, so it is better placed at the base.



Axis 500 – 700 Premium Line

Mobile Operating Tables with most
cantilevered column for all surgical
applications with dedicated Accessoires



Axis – Technical Specifications

Axis - Technical Specifications	Axis 500	Axis 600	Axis 700
Table dimensions incl. head- and leg plates with side rails (mm)	L 1997 x W 580	L 1997 x W 580	L 485 x W 580
X-Ray window (mm)	1450	1450	1450
Height adjustment with pad (mm)	755 - 1161	755 - 1161	755 - 1161
Height adjustment without pad (mm)	686 - 1092	686 - 1092	686 - 1092
Head-down Trendelenburg position	28°	28°	28°
Leg-down/Reverse-Trendelenburg position	28°	28°	28°
Seat plate up (motorized)	76°	76°	76°
Seat plate down (motorized)	40°	40°	40°
Leg plate up (manual) (Arcus 601 motorized)*	15°	15°+(15°electric)	Optional 15°
Leg plate down (manual) (Arcus 601 motorized)*	95°	95°+(70°electric)	Optional 95°
Head plate up (manual)	30°	30°	30°
Head plate down (manual)	30°	30°	30°
Lateral tilting to both sides	18°	18°	18°
Side rail dimensions (mm)	25 x 10	25 x 10	25 x 10
Dimension X-ray cassette channel (mm)	W 470 x H 33	W 470 x H 33	W 470 x H 33
Dimension X-ray tray/max. size cassetts (mm)	L 315 x W 403	L 315 x W 403	L 315 x W 403
Table weight without head and leg plates (kg)	284	292	271
Max. patient weight capacity (kg)	350	350	350
Dynamic weight capacity, all positions (kg)	220	220	220
Pad thickness (mm)	69	69	69
Wheels dimensions (mm)	120	120	120
IPX – Protection against entry of liquids	IPX4	IPX4	IPX4




Segment Movements	Axis 500	Axis 600	Axis 700
Height adjustment up/down	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Trendelenburg / Reverse trendelenburg	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Lateral tilt left / right	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Seat plate up / down	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Leg plate up / down	Manual	Electro-hydraulic In A - position	Manual
Head plate up / down	Manual	Electro-hydraulic In B - position	Manual
Flex / Reflex	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Zero position	Electro-hydraulic	Electro-hydraulic	Electro-hydraulic
Leg plates swivelling out 90°	Manual	Manual	Manual

Electrical Data


- Electrical mains connection 100-120 / 200-240 V AC, 50-60 HZ
- Control voltage 24 V DC
- Main system: (2 x 12 V DC) capacity 12 AH
- Emergency system: 24 V DC (2 x 12 V DC) capacity 7 AH
- 4 maintenance-free batteries
- Battery charging time approx. 6 to 8 hours
- Battery capacity approx. 20 operations
- Charge indication on the hand control
- Low battery acoustic signal




Axis – Trendelenburg / Reverse Trendelenburg

Degrees	Highest Position	Lowest Position
		


Axis – Head Rest




Axis – Leg Plates




Axis – Beach-Chair-Position




Axis – Flex




Axis – Reflex




Axis – Lateral Tilt



Axis – Dimensions



Axis 600 – With Second Motorized Joint



Axis – Technical Specifications

Axis 500 - Segments



Axis 600 - Segments



Axis - Range



Without pads: 686 to 1092 mm

With pads: 755 to 1161 mm

Height pads: 69 mm


Axis - Patient Weight Capacity

350 kg
Static


220 kg
Dynamic



Axis – X-Ray Window



Axis - Foot Dimensions



Axis - Controls

Hand Control



This Hand control with cable allows to realize all the movements of the table, and in addition it has a function that indicates all the errors that can occur, by system of codes that allows a fast and effective detection.

Emergency Control




Emergency control that allows all movements to be carried out with a second system that activates the hydraulic system, in the event that the control stops working.









Foot Switch






Item:90424
Foot switch for **Axis 500-700**

Axis – Components and Carbon Plates

Axis 500		 ➤ Removable seat plate ➤ Can be used for all surgeries ➤ Can attach all accessories CSF / MTS or any compatible Schaerer traction device
Item: 91595 Axis 500		
Item: 91260 Packaging		
Axis 600		 ➤ Seat plate not removable ➤ Ideally for Gyn/Obs, Urology and Neurosurgery ➤ Not for Ortho /Trauma
Item: 91596 Axis 600		
Item: 91260 Packaging		
Axis 700		 ➤ Table same like Axis 500 but without seat plate ➤ Can be used for all surgeries ➤ Can attach all accessories CSF / MTS or any compatible Schaerer traction device
Item: 91597 Axis 700		
Item: 91260 Packaging		



Head Plates		Table Extension
Item: 91490 Head rest standard 	Item: 91491 Head rest rigid 	Item: 90939 Table extension for Axis 
Leg Plates		
Item: 91495 Leg plates 2-part 	Item: 91496 Leg plates 2-part triangular 	Item: 90914 Leg plate one-piece 
Intermediate Plates		X-Ray Holder
Item: 90952 Intermediate plate with Gyn cut-out 	Item: 90951 Intermediate plate 	Item: 90908 X-ray holder for Axis 500 Item: 90909 X-ray holder for Axis 600 
Carbon Plates		
Item: 90627 Carbon plate 1400mm for Axis 500/700 (Instead of Seat plate) 	Item: 90624 Carbon plate 850 mm for Axis 500/700 (Instead of Seat plate) 	Item: 90625 Carbon plate 850 mm for Axis 500/600 (Instead of Leg plates) 

Combinations	Axis – Components Connection Possibilities					Length			
Combination A	A Position	Head Plate 	+	Axis 500/600 	+	Leg Plates 	= 1997 mm		
Combination A	B Position	Leg Plates 	+	Axis 500/600 	+	Head Plate 	= 1997 mm		
Combination B		Head Plate 	+	Int. Plate Item: 90951 	+	Axis 500/600 	+	Leg Plates 	= 2370 mm
Combination C		Head Plate 	+	Int. Plate Item: 90952 	+	Axis 500/600 	+	Leg Plates 	= 2370 mm
Combination D		Head Plate 	+	Axis 500/600 	+	Int. Plate Item 90951 		Leg Plates 	= 1645 mm
Combination E		Head Plate 	+	Axis 500/600 	+	Int. Plate Item 90952 		Leg Plates 	= 1645 mm
	Axis – Carbon Plates Connection Possibilities								
Combination F	A - Position	Head Plate 	+	Axis 500/600 	+	Item:90625 			= 2122 mm
Combination G	A - Position	Head Plate 	+	Axis 500/700 	+	Item:90627 			= 2284 mm
Combination H		Head Plate 	+	Int. Plate Item: 90951 	+	Axis 500/700 	+	Item:90627 	= 2657 mm
Combination I		Head Plate 	+	Axis 500/700 	+	Item:90624 			= 1734 mm
Combination J		Head Plate 	+	Int. Plate Item: 90951 	+	Axis 500/700 	+	Item:90624 	= 2107 mm
Combination K		Head Plate 	+	Int. Plate Item: 90952 	+	Axis 500/700 	+	Item:90624 	= 2107 mm

Accessories for General Surgery

6  Item: 90321 (or)  Item: 90323 (+)  Item: 90506 (or)  Item: 90310  Item: 90635 (or)  Item: 90636 (+)  Item: 90506 (or)  Item: 90310	7  Item: 90520 (or)  Item: 4646.05 (+)  Item: 5.890-VE 8  Item: 800-0057 9  Item: 92021 10  Item: 90500	11  Item: 90332 12  Item: 90360 13  Item: 92028 14  Item: 92027 15  Item: 90913 16  Item: 90528	<table><tr><th>No</th><th>Item</th><th>Qty</th><th>Description</th></tr><tr><td rowspan="2">6</td><td>90321</td><td>1</td><td>Anaesthesia screen, rigid (+)</td></tr><tr><td>90506</td><td>1</td><td>Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp</td></tr><tr><td rowspan="2">6</td><td>90323</td><td>1</td><td>Anaesthesia screen with central articulation (+)</td></tr><tr><td>90506</td><td>1</td><td>Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp</td></tr><tr><td rowspan="2">6</td><td>90635</td><td>1</td><td>Anaesthesia screen, with extendable rods (+)</td></tr><tr><td>90506</td><td>1</td><td>Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp</td></tr><tr><td rowspan="2">6</td><td>90636</td><td>1</td><td>Anaesthesia screen central art., ext. rods (+)</td></tr><tr><td>90506</td><td>1</td><td>Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp</td></tr><tr><td rowspan="3">7</td><td rowspan="3">90520</td><td>2</td><td>Arm rest (or)</td></tr><tr><td></td><td>4646.05 Arm support quick central fixation (+)</td></tr><tr><td></td><td>5.890-VE Fastening straps (1 pair)</td></tr><tr><td>8</td><td>800-0057</td><td>1</td><td>Restraint Strap</td></tr><tr><td>9</td><td>92021</td><td>1</td><td>Arm strap with fixing screw</td></tr><tr><td>10</td><td>90500</td><td>1</td><td>Infusion pole with clamp, height-adjustable</td></tr><tr><td>11</td><td>90332</td><td>1</td><td>Instrument tray with clamp, height-adjustable</td></tr><tr><td rowspan="2">12</td><td>90360</td><td>1</td><td>Accessories trolley (+)</td></tr><tr><td>91300</td><td>1</td><td>Packaging for Accessories trolley</td></tr></table> <table><tr><th colspan="4">Dedicated Accessories</th></tr><tr><td>13</td><td>92028</td><td>1</td><td>Universal Adapter for DORO® basic unit</td></tr><tr><td>14</td><td>92027</td><td>1</td><td>Collecting tray TUR for Axis 500/600</td></tr><tr><td>15</td><td>90913</td><td>1</td><td>Hand surgery table with fixing screw</td></tr><tr><td>16</td><td>90528</td><td>1</td><td>X-ray cassette frame for Axis 500/600</td></tr></table>	No	Item	Qty	Description	6	90321	1	Anaesthesia screen, rigid (+)	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp	6	90323	1	Anaesthesia screen with central articulation (+)	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp	6	90635	1	Anaesthesia screen, with extendable rods (+)	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp	6	90636	1	Anaesthesia screen central art., ext. rods (+)	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp	7	90520	2	Arm rest (or)		4646.05 Arm support quick central fixation (+)		5.890-VE Fastening straps (1 pair)	8	800-0057	1	Restraint Strap	9	92021	1	Arm strap with fixing screw	10	90500	1	Infusion pole with clamp, height-adjustable	11	90332	1	Instrument tray with clamp, height-adjustable	12	90360	1	Accessories trolley (+)	91300	1	Packaging for Accessories trolley	Dedicated Accessories				13	92028	1	Universal Adapter for DORO® basic unit	14	92027	1	Collecting tray TUR for Axis 500/600	15	90913	1	Hand surgery table with fixing screw	16	90528	1	X-ray cassette frame for Axis 500/600
No	Item	Qty	Description																																																																																			
6	90321	1	Anaesthesia screen, rigid (+)																																																																																			
	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp																																																																																			
6	90323	1	Anaesthesia screen with central articulation (+)																																																																																			
	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp																																																																																			
6	90635	1	Anaesthesia screen, with extendable rods (+)																																																																																			
	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp																																																																																			
6	90636	1	Anaesthesia screen central art., ext. rods (+)																																																																																			
	90506	1	Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp																																																																																			
7	90520	2	Arm rest (or)																																																																																			
			4646.05 Arm support quick central fixation (+)																																																																																			
			5.890-VE Fastening straps (1 pair)																																																																																			
8	800-0057	1	Restraint Strap																																																																																			
9	92021	1	Arm strap with fixing screw																																																																																			
10	90500	1	Infusion pole with clamp, height-adjustable																																																																																			
11	90332	1	Instrument tray with clamp, height-adjustable																																																																																			
12	90360	1	Accessories trolley (+)																																																																																			
	91300	1	Packaging for Accessories trolley																																																																																			
Dedicated Accessories																																																																																						
13	92028	1	Universal Adapter for DORO® basic unit																																																																																			
14	92027	1	Collecting tray TUR for Axis 500/600																																																																																			
15	90913	1	Hand surgery table with fixing screw																																																																																			
16	90528	1	X-ray cassette frame for Axis 500/600																																																																																			

Ophthalmology / ENT – Configuration A



Ophthalmology / ENT – Configuration B



No	Item	Qty	Description
1	92028	1	Adapter for DORO® basic unit, for Axis
2	90045	1	Horseshoe-Head rest holder 3 joints
3	1002.72A0	1	Horseshoe-shaped Head rest, two pieces
4	1004.75A0	1	Wrist hand support for Ophthalmology

Abdominal – Configuration A

1



(or) (+)

Item: 90520 Item: 4646.05
Item: 5.890-VE

2



(+) Item: 90506 (or)
Item: 90310

3



(+) Item: 90506 (or)
Item: 90310



No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90325 90506	1 1	Arm strap with fixing screw Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	90330 90506	2 2	Leg holder "Goepel" (1 piece) (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp

Abdominal – Configuration B

1



(or) (+)

Item: 90520 Item: 4646.05
Item: 5.890-VE

2




(+) Item: 90506 (or)
Item: 90310

3




(+) Adjustable Leg holder



No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90325 90506	1 1	Arm strap with fixing screw Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	800-0342-PR 800-0338-EU	1 2	Adjustable Leg holder for patients up to 182 kg (or) 800-0342-PL for patients up to 272 kg (or) 800-0364 for patients up to 454 kg (+) Springloc Clamp (order 2 pieces)

Bariatric, Laparoscopic – Configuration A

1  Item: 90520 (or)  Item: 4646.05 + 5.890-VE	2  Item: 90079 (+)  Item: 81501	3  Item: 800-0055 (+)  Item: 800-0338-EU	4  Item: 90939
--	--	--	---



No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90079 81501	1 2	Knee cradles (1 pair), with pads (+) Clamp, Knee cradles (order 2 pieces)
Optional			
3	800-0055 800-0338-EU	1 2	Nissen straps, 272 kg patient w. cap. (1 pair) (+) Springloc Clamp
4	90939	4	Table extension for Axis (1 piece),

Bariatric, Laparoscopic – Configuration B

1  Item: 90520 (or)  Item: 4646.05 + 5.890-VE	2  (+) Adjustable Leg holder	3  Item: 800-0055 (+)  Item: 800-0338-EU	4  Item: 90939
--	--	--	---



No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	800-0342-PR 800-0338-EU	1 2	Adjustable Leg holder for patients up to 182 kg (or) 800-0342-PL for patients up to 272 kg (or) 800-0364 for patients up to 454 kg (+) Springloc Clamp (order 2 pieces)
3	800-0055 800-0338-EU	1 2	Nissen straps, 272 kg patient w. cap. (1 pair) (+) Springloc Clamp
4	90939	4	Table extension for Axis (1 piece),

Nephrectomy, Thorax

1



Item: **90520**

(or)



Item: **4646.05**
Item: **5.890-VE**

(+)



2



Item: **90325**

(+)



Item: **90506**

(or)



Item: **90310**

3



Item: **90351**

(+)



Item: **90350**

(+)



Item: **90539**

(or)



Item: **90314**

4



Item: **90352**

(+)



Item: **90350**

(+)



Item: **90539**


(or)



Item: **90314**

2

3



1

5



Item: **90609**

(+)



Item: **90350**

(+)




Item: **90539**

(or)



Item: **90314**



No	Item	Qty	Description
1	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90325 90506	1 1	Arm rest high, without clamp (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	90351 90350 90539	1 1 1	Side and back rest, with pad (+) Universal holder for side and back rest (+) Radial clamp with fixing screw, rotating (or) 90314 Rapid-action radial setting clamp
4	90352 90350 90539	1 1 1	Side and back supporting roll Universal holder for side and back rest (+) Radial clamp with fixing screw, rotating (or) 90314 Rapid-action radial setting clamp
5	90609 90350 90539	1 1 1	Sacrum/Pubis/Sternum support, with pad (+) Universal holder for side and back rest (+) Radial clamp with fixing screw, rotating (or) 90314 Rapid-action radial setting clamp

Rectal, Andrews, Genupectoral

1



Item: 90520

Item: 4646.05
Item: 5.890-VE

2




Item: 90277

Item: 90506

Item: 90310


3



Item: 90079



Item: 81501

4



Item: AP022

Item: QPHPD-ADULT 1



No	Item	Qty	Description
1	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90277 90506	1 2	Abdominal rest, with pad, adjustable (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	90079 81501	1 2	Knee cradles (1 pair), with pads (+) Clamp, Knee cradles (order 2 pieces)
4	AP022	1	Gel Headrest for prone position (or) QPHPD-ADULT 1 – Head rest (single use)


Gynaecology / Urology – Configuration A



Item: 90330



Item: 90421



Item: 90520
(or)
Item: 4646.05
Item: 5.890-VE



No	Item	Qty	Description
1	90330 90506	2 2	Leg holder "Goepel" (1 piece) (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
2	90421	1	Collecting tray TUR for Axis (or)
3	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)

Gynaecology / Urology – Configuration B



Adjustable Leg holder



Item: 90421



Item: 90520
(or)
Item: 4646.05 + 5.890-VE



Item: 800-0074-S



No	Item	Qty	Description
1	800-0342-PR 800-0342-PL 800-0364 800-0338-EU	1 1 1 2	Adjustable Leg holder for patients up to 182 kg (or) Adjustable Leg holder for patients up to 272 kg (or) Adjustable Leg holder for patients up to 454 kg (+) Springloc Clamp (order 2 pieces)
2	90421	1	Collecting tray TUR for Axis
3	90520	2	Arm rest with clamp (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
Optional			
4	800-0074-S	1	Great White® Stirrup Dolly

CFP Carbon Fibre Plates

Item: 90627


Carbon-fibre plate 1400 mm for
Axis 500/700

(Instead of seat plate)

Patient weight capacity: 130kg

Dimensions: 1400 x 520 x 130 mm

(not for Axis 600)



Configuration

No	Item	Qty	Description
1	91597 91595	1 pc	Axis 700 (or) Axis 500 (Without seat plate)
2	91490	1 pc	Head rest standard
3	90627	1 pc	Carbon-fibre plate 1400 mm
4	90593	2 pc	Side rail 25 x 10 mm (1 piece)
5	AP004	1 pc	Head support, 20 x 7.5 x 4.5 cm.
6	91260	1 pc	Packaging

Item: 90624


Carbon-fibre plate 850 mm for
Axis 500/700

(Instead of seat plate)

Patient weight capacity: 220kg

Dimensions: 850 x 520 x 130 mm

(not for Axis 600)



Configuration

No	Item	Qty	Description
1	91597 91595	1 pc	Axis 700 (or) Axis 500 (Without seat plate)
2	91490	1 pc	Head rest standard
3	90624	1 pc	Carbon-fibre plate 850 mm
4	90593	2 pc	Side rail 25 x 10 mm (1 piece)
5	AP004	1 pc	Head support, 20 x 7.5 x 4.5 cm.
6	91260	1 pc	Packaging

Item: 90625


Carbon-fibre plate 850 mm for
Axis 500/600

(Instead of the leg plate)

Patient weight capacity: 90kg

Dimensions: 850 x 520 x 130 mm

(not for Axis 700)



Configuration



No	Item	Qty	Description
1	91595 91596	1 pc	Axis 500 (or) Axis 600
2	91490	1 pc	Head rest standard
3	90625	1 pc	Carbon-fibre plate 850 mm
4	90593	2 pc	Side rail 25 x 10 mm (1 piece)
5	AP004	1 pc	Head support, 20 x 7.5 x 4.5 cm.
6	91260	1 pc	Packaging




For example: Axis 700 + Head rest + Carbon plate (1400 mm) = 2284 mm



Spine Surgery – Configuration

1  Item: NO 165	2  Item: NA 401	3  Item: 90520 Item: 4646.05 Item: 5.890-VE
4  Item: 92028  Item: 90456  Item: 90041  Item: 90003 (or)  Item: QPHPD-ADULT 1 (or)  Item: AP022		
5  Item: 90625	6  Item: 800-0235	7  Item: 800-0074-SF



No	Item	Qty	Description
1	NO 165	1	Torso Pad "Onsernone"
2	NA 401	1	Leg support wave cushion 65x50x13cm
3	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
4	92028 90456 90041 90003	1 1 1 1	Adapter for DORO® basic unit, for Axis (+) Hexagonal adapter with square adapter (+) Head rest holder, with lateral adjustment (+) Head rest (Horse shoe) with pad (or) QPHPD-ADULT 1 - QView™ Head positioning device (single use) (or) AP022 Head rest for prone position
Optional Accessories			
5	90625	1	Carbon-fibre plate 850 mm
6	800-0235	1	Schure Spine Frame, Radiolucent
7	800-0074-SF	1	Schure Spine Frame Dolly

CSF Carbon Spine Frame

1



Item: 91600

2



Item: 90520


Item: 4646.05
Item: 5.890-VE

4



Item: 91603

(or)



Item: 91604

5



Item: AP022

(or)



Item: QPHPD-ADULT 1

Item: 91600


Schaerer Carbon Spine Frame

Patient weight capacity:

Axis: 200 kg



No floor
support

3



360°
Radiolucent


No	Item	Qty	Description
1	91600	1	Schaerer Carbon Spine Frame (CSF) for Axis 500/700
2	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
3	90776	1	Stability support with lock. lever for Axis, removable
Optional Accessories			
4	91603 91604	1	Side rail adapter short (or) Side rail adapter long
5	AP022	1	Head rest for prone position (or) QPHPD-ADULT 1 - QViewTM Head positioning device (single use)

CSF - Carbon Spine Frame DORO® Radiolucent Head Rest System






(or)



Radiolucent			
No	Item	Qty	Description
6	1101.021	1	DORO® LUCENT Basic support
7	1101.031	1	DORO® LUCENT Clamp intermediate piece
8	1101.026	1	DORO® LUCENT Intermediate piece
9	3033-00	1	DORO® Swivel adapter, radiolucent
10	3034-00	1	DORO® Skull clamp, radiolucent Quick-Rail



Aluminium			
No	Item	Qty	Description
11	3001-00	1	DORO® Base unit aluminum
12	3002-00	1	DORO® Swivel adaptor aluminum (or)
13	1001.001	1	DORO® Skull clamp QR3 aluminum

Technical Data




Patient weight capacity	200 kg
Overall length	1205 mm
Length 360° radiolucent	805 mm
Frame length	1052 mm
Width frame	480 mm

Neuro – Configuration A - Aluminium Skull Clamp System



No	Item	Qty	Description
1	92028	1	Adapter for DORO® Basic unit, Axis (or) 3011-00 Rail adapter (for all tables)
2	3001-00	1	DORO® Base unit aluminum
3	3002-00	1	DORO® Swivel adaptor aluminum (or) DORO® Ball adapter, allows 360° rotation
4	1001.001	1	DORO® Skull clamp QR3 aluminum
5	3009-00	1	DORO® Horse shoe head rest for adults
6	3008-00	1	DORO® Horse shoe head rest for children
		1	Pins reusable (or) Pins disposable
Additional accessories for sitting position			
7	3007-00 90506	1	DORO® Rail adapter for sitting position (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp



Pins reusable



Item: 3005-00 (Stainless Steel)
Item: 3005-50 (Titanium)
Reusable skull pins for **adults**
(3 pieces)



Item: 3004-00 (Stainless Steel)
Item: 3004-50 (Titanium)
Reusable skull pins for **children**
(3 pieces)

Pins disposable



Item: 1106.003
DORO® disposable skull pins, for **adults**,
(36 pieces)



Item: 3006-10
DORO® disposable skull pins, for **children**,
(36 pieces)

Pins disposable



Item: 3006-20
DORO® disposable skull pins, for **adults**,
(36 pieces)



Item: 3006-30
DORO® disposable skull pins, for **children**,
(36 pieces)

Neuro – Configuration B - Radiolucent Skull Clamp System

Item: **1101.020** (Complete Set)

DORO® Radiolucent Head Rest System Compatible with all tables
(No adapter needed)



Item: **1101.021**



Item: **1101.031**




Item: **1101.026**



Item: **3033-00**





Item: **3034-00**



No	Item	Qty	Description
1	1101.021	1	DORO® LUCENT Basic support
2	1101.031	1	DORO® LUCENT Clamp intermediate piece
3	1101.026	1	DORO® LUCENT Intermediate piece
4	3033-00	1	DORO® Swivel adapter, radiolucent
5	3034-00	1	DORO® Skull clamp, radiolucent
6		1	Pins reusable (or) Pins disposable
Additional accessories for sitting position			
7	3007-00 90506	1	DORO® Rail adapter for sitting position (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
8	3001-00	1	DORO® Base unit aluminum

Sitting Position



Pins reusable



Item: **3005-00** (Stainless Steel)
Item: **3005-50** (Titanium)
Reusable skull pins for **adults**
(3 pieces)



Item: **3004-00** (Stainless Steel)
Item: **3004-50** (Titanium)
Reusable skull pins for **children**
(3 pieces)

Pins disposable



Item: **1106.003**
DORO® disposable skull
pins, for **adults**,
(36 pieces)



Item: **3006-10**
DORO® disposable
skull pins, for **children**,
(36 pieces)

Pins disposable



Item: **3006-20**
DORO® disposable
skull pins, for **adults**,
(36 pieces)



Item: **3006-30**
DORO® disposable
skull pins, for **children**,
(36 pieces)

MTS Modular Traction System (One Surgical Table - Multiple Disciplines)


Applications

<u>General Surgery</u>	<u>Femur</u>	<u>Tibia</u>	<u>Hip</u>	<u>Lateral Position</u>	<u>Shoulder</u>
					


Dimensions in retracted position (in mm)

Dimensions in extended position (in mm)


- MTS




- MTS + MIS



- MTS












- MTS + MIS





Dimensions seat plate (in mm)


Dimensions MTS-Trolley (in mm)






MTS - Femur			
No	Images	Item	Description
1		90962	Seat Plate ADVANCED, Axis 500/700
2		90989	Supplement bar with plate, left, Axis 500/700
3		90988	Supplement bar with plate, right, Axis 500/700
4		90971	Perineal rod
5		90990	Perineal rod for arthroscopy
6		90958	Offset holder
7		90983	Additional joint
8		90981	Traction bar Metal
9		90974	Traction bar ADVANCED

MTS - Femur			
No	Images	Item	Description
10		90978	Multiclamp (For Traction bar ADVANCED only)
11		90928	Multiclamp slide (For Traction bar ADVANCED only)
12		90980	Traction unit
13		90895	Traction unit, vertical shift
14		90763	Extension shoe MTS, white (1 piece)
15		RO.0162 .2018	Traction shoe with side rail interface
16		RO.0148 .2018	Pad set incl. fixation straps (10 pieces)
17		90330 (+) 90506	Leg holder "Goepel" (1 piece) (+) 90506 Radial setting clamp with fixing screw
18		90994	Extension bow holder, for Kirschner- wire or pin

MTS – Tibia additional Items

No	Images	Item	Description
19		90780	Central rail adapter
20		90897	Radial setting clamp with star screw
21		90759	Holder for knee rest
22		90896	Adapter for lower leg surgery
23		92699 92698	Foot plate Adult Foot plate for child
24		90193	Telescopic support



MTS – Shoulder additional Items

No	Images	Item	Description
25		90762	Head rest system, Axis 500/700
26		F-90023	Head Pad for MTS (10 pieces)
27		90948	Extension plate, Axis 500/700
28		90925	Lateral ext. plate for head rest system, left, Axis 500/700
29		90924	Lateral ext. plate for head rest system, right, Axis 500/700

MTS – Lateral position additional Items

No	Images	Item	Description
30		92611	ETS Counter Traction Unit for lateral position

Optional Items

31		90931	Gyn. pad set
32		90959	Gyn. cut-out plate, Axis 500/700





MTS - MIS Items

No	Images	Item	Description
34		90738 91298	Bar MIS for MTS system (1 piece) Packaging MIS
35		90775 91299	MIS Trolley Packaging Trolley

MTS - Trolley




No	Images	Item	Description
33		90957 91304	MTS Trolley Packaging Trolley Measurements: Width: 625 mm Length: 890 mm Height: 1610 mm

MTS – Table Configuration




No	Item	Qty	Description
A	91597	1	Axis 700
B	91490	1	Head rest standard
C	91495 91496	1	Leg plates 2-part (or) Leg plates 2-part (triangular)
36	90776	1	Stability support
Basic MTS Accessories			
1	90962	1	Seat Plate ADVANCED, Axis 500/700
2	90989	1	Supplement bar with plate, left, Axis 500/700
3	90988	1	Supplement bar with plate right, Axis500/700

Femur - Traction with Metal bars (Additional to MTS Table)


4  Item: 90971	6  Item: 90958	7  Item: 90983	8  Item: 90981	
12  Item: 90980	(or) 13  Item: 90895	14  Item: 90763	(or) 15  Item: RO.0162.2018 RO.0148.2018	16 

1 Leg Configuration



No	Item	Qty	Description
4	90971	1	Perineal rod
6	90958	1	Offset holder (Optional)
7	90983	1	Additional joint (Optional)
8	90981	1	Traction bar Metal
12	90980	1	Traction unit (or)
13	90895	1	Traction unit, vertical shift
14	90763	1	Extension shoe MTS, white (1 piece) (or)
15	RO.0162.2018	1	Traction shoe with side rail interface (+)
16	RO.0148.2018	1	Pad set straps (disposable) (10 pieces)
17	90330 90506	1	Leg holder "Goepel" (1 piece) (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp


2 Leg Configuration



No	Item	Qty	Description
4	90971	1	Perineal rod
6	90958	1	Offset holder (Optional)
7	90983	1	Additional joint (Optional)
8	90981	2	Traction bar Metal
12	90980	2	Traction unit (or)
13	90895	2	Traction unit, vertical shift
14	90763	2	Extension shoe MTS, white (1 piece) (or)
15	RO.0162.2018	2	Traction shoe with side rail interface (+)
16	RO.0148.2018	2	Pad set straps (disposable) (10 pieces)


Femur - Traction with Carbon bars (Additional to MTS Table)

1 Leg Configuration




No	Item	Qty	Description
4	90971	1	Perineal rod
6	90958	1	Offset holder (Optional)
7	90983	1	Additional joint (Optional)
9	90974	1	Traction bar ADVANCED
10	90978	1	Multiclamp (or)
11	90928		Multiclamp slide
12	90980	1	Traction unit (or)
13	90895		Traction unit, vertical shift
14	90763	1	Extension shoe MTS, white (1 piece) (or)
15	RO.0162.2018		Traction shoe with side rail interface (+)
16	RO.0148.2018		Pad set straps (dispos.) (10 pieces)
17	90330	1	Leg holder "Goepel" (1 piece) (+)
	90506		Radial setting clamp with screw (or)
	90310		Rapid-action radial setting clamp

2 Leg Configuration



No	Item	Qty	Description
4	90971	1	Perineal rod
6	90958	1	Offset holder (Optional)
7	90983	1	Additional joint (Optional)
9	90974	2	Traction bar ADVANCED
10	90978	2	Multiclamp (or)
11	90928		Multiclamp slide
12	90980	2	Traction unit (or)
13	90895		Traction unit, vertical shift
14	90763	2	Extension shoe MTS, white (1 piece) (or)
15	RO.0162.2018		Traction shoe with side rail interface (+)
16	RO.0148.2018		Pad set straps (dispos.) (10 pieces)

Tibia - in addition to Femur Traction with **Carbon** or **Metal** bars

19	20	21	22	23	24
					
Item: 90780	Item: 90897	Item: 90759	Item: 90896	Item: 90699	Item: 90193



No	Item	Qty	Description
19	90780	1	Central rail adapter
20	90897	1	Radial setting clamp with star screw
21	90759	1	Holder for knee rest
22	90896	1	Adapter for lower leg surgery
23	92699	1	Foot plate / 92698 foot plate for child
24	90193	1	Telescopic support

Lateral Position - in addition to Femur Traction with **Carbon** or **Metal** bars




Item: **92611**



Item: **90325 (+) 90506**







Item: **90520**



No	Item	Qty	Description
30	92611	1	Counter Traction Unit for lateral position
	90325	1	Arm rest high, without clamp (+) 90506
	90520	1	Arm rest (or) (4646.05 (+) 5.890-VE)

Shoulder (Additional to MTS Table)

25  Item: 90762	26  Item: F-90023	27  Item: 90948	28  Item: 90925	29  Item: 90924
---	---	---	--	---



No	Item	Qty	Description
25	90762	1	Head rest system, Axis 500/700
26	F-90023	1	Head Pad for MTS (10 pieces)
27	90948	1	Extension plate, Axis 500/700
28	90925	1	Lateral ext. plate for head rest system, left, Axis 500/700
29	90924	1	Lateral ext. plate for head rest system, right, Axis 500/700

Possible to combine Metal with Carbon

MTS-MIS



Item: **90738**
Bar MIS, for MTS



Hand Surgery



1



1

(or)



2

(+)



Item: **90506**
(order 2 pc.)

No	Item	Qty	Description
1	90913	1	Hand surgery table, with fixing screw (or)
2	800-0040	1	Hand surgery table, radiolucent
	90506	2	Radial setting clamp

Humerus / Traction Device



1





No	Item	Qty	Description
1	9907038	1	Humerus and Forearm Positioning Device

Knee Arthroscopy



1

(or)



2

(+)



Item: **90539**



Item: **90314**



2

(+)



Item: **90506**



Item: **90310**





1



2

No	Item	Qty	Description
1	90357	1	Leg holder arthroscopy (or)
	90356	1	Knee rest (+)
	90539	1	Radial clamp with fixing screw, rotating (or)
	90314	1	Rapid-action radial clamp, rotating
2	90330	1	Leg holder "Goepel" (1 piece) (+)
	90506	1	Radial setting clamp with fixing screw (or)
	90310	1	Rapid-action radial setting clamp



SSP - Schaerer Shoulder Plate - Configuration A


1  Item: 92645	2  Item: 90456	3  Item: 90041	4  Item: 90566 (or) Item: 90003
--	--	---	---



No	Item	Qty	Description
1	92645	1	Shoulder plate (SSP) for Axis 500-700
2	90456	1	Adapter
3	90041	1	Head rest holder, with lateral adjustment
4	90566 90003	1	Helmet, radiolucent (or) Head rest (Horse shoe) with pad

SSP - Schaerer Shoulder Plate - Configuration B

1  Item: 92645	2  Item: 92666	3  Item: 90899	4  Item: F-90023
--	--	---	--



No	Item	Qty	Description
1	92645	1	Shoulder plate (SSP) for Axis 500-700
2	92666	1	Head rest holder
3	90899	1	Head rest
4	F-90023	1	Head pad (10 Pieces) (single use)

1 - Highest capacity for eccentric column table

Highest weight capacity for an eccentric column table. (350kg/220kg) means more safety for all kind of patient positions, all kind of operations, all surgical disciplines and a much more stable table. There is no sliding who limits the weight capacity



2 - Highest-grade stainless-steel table

Whole table is made of highest-grade stainless steel. Special treatment of all stainless-steel surfaces to withstand the most aggressive disinfection liquids. Table fully covered; no open spaces left.



3 - Largest X ray window in the market

Due to its very eccentric column design, highest X ray window in market. No need for sliding! From others copied but never achieved!



4 - Highest Durability of hand control

Hand control design with several functions and LCD screen; Operating functions of the table display describes movements of the table that allows to control movements even when the table is fully draped and can't be seen. Automatic system check up show failures for a 1st level service info without delay. Programming features, different languages, adjust "Level" position, additional check up of connections.




5 - Highest safety features

Axis 500-700 is additionally equipped with a second built-in panel controlling an independent circuit. The (override) panel controls a second pump and motor which is simultaneous online and can be used instead of the main system at any time, even if there is no failure of the main system or hand control. If in an emergency situation, the main system fails, the operating table can be driven with the emergency panel.




6 - Highest patient comfort

Our viscoelastic pads are electric conductive, silicon and phthalate free, 69 mm high, Velcro free, ultrasound welded, without stitches. easily removable and attachable. Appropriate for long operations and heavy weight patients.



7 - Highest mobility and steerability

High grade durable double swiveling castors, for easy steering and maneuvering four electrohydraulic driven floor locks secure the table in its position . Lock and release by the electric foot locking lever. Steering is easy and convenient by the Directional lock



AXIS 350 **Advanced Line**

Mobile Operating Table with sliding top
for all surgical applications with wide
range of Accessories



Axis 350 - Specifications	Specifications
Length with Head Rest, Long Intermediate Plate and Leg Plates	2185 mm
Length with Head Rest, Short Intermediate Plate and Leg Plates	2035 mm
Width of table top (without side rails)	520 mm
Width of table top (with side rails)	580 mm
Length table pad (without dedicated accessories)	630 mm
Width table pad (without dedicated accessories)	520 mm
Column moving range (high & low)	370 mm
Sliding table top	300 mm
Column Range axis 350 High (lowest & highest position, without pad)	750 - 1120 mm
Column Range axis 350 Low (lowest & highest position, without pad)	612 - 982 mm
Pad thickness	69 mm
X-ray window A-Position	1285 x 348 mm
X-ray window B-Position	1135 x 348 mm
Side rail dimension	25 x 10 mm
Table weight (without dedicated accessories)	255 kg
Static workload	450 kg
Dynamic workload	220 kg
Wheels diameter	125 mm
Seat Plate up / down	80° / 40°
Trendelenburg / Reversed Trendelenburg	30°
Lateral tilting left / right	20°
Head Rest up / down manual	45°
Leg Plate up / down manual	-90° / 20°
IPX – Protection against entry of liquids	IPX4

Batteries

➤ Main Power Supply

- ✓ 2 Rechargeable Lead Acid Batteries
- ✓ 22 Ah, 12V (24V as System)
- ✓ Lifetime: 500 Cycles approx. 2 Years
- ✓ 50 Operations, capacity of appr. One week usage

➤ Emergency Power Supply

- ✓ 2 Rechargeable Lead Acid Batteries
- ✓ 3.2 Ah, 12V (24V as System)
- ✓ Lifetime: 500 Cycles approx. 2 Years

➤ MOP Power Supply

- ✓ Rechargeable Lithium-Ionen Battery
- ✓ 2.6 Ah, 3.6V
- ✓ Lifetime: 500 Cycles approx. 2 Years




❖ For long battery life – charge table every day over night!

Axis 350 - Dimensions

High and Low Dimension	Sliding	Head plate
Trendelenburg / Reverse Trendelenburg	Lateral tilt left / right	Leg plate (One-piece & 2-part Leg Plates)
Flex / Reflex-Position		Foot dimensions


Components - Base

- 1 - Head plate
- 2 - Intermediate plate
- 3 - Leg plate
- 4 - Emergency operational panel
- 5 - Socket for wired hand control
- 6 - Lever for foot pump
- 7 - Terminal for potential equalization
- 8 - Socket for control device
- 9 - Socket for main cable
- 10 - Battery charging indicator



Axis 350 – Technical Specifications

Axis 350 – Segments



Axis 350 – Sliding Top

Manual sliding table top **300 mm**
with **electric** lock/unlock button



Axis 350 - Range

Low Version

- Without pads: 612 to 982 mm
- With pads: 681 to 1051 mm
- Height pads: 69 mm



High Version

- Without pads: 750 to 1120 mm
- With pads: 819 to 1189 mm
- Height pads: 69 mm



Axis 350 – Length

Head Rest, **Long** Intermediate Plate and Leg Plates



Head Rest, **Short** Intermediate Plate and Leg Plates



Axis 350 – Patient Weight Capacity

Static

450 kg
Static



Dynamic **A-Position** full slide

185 kg
Dynamic



Dynamic **B-Position** full slide

220 kg
Dynamic



Axis 350 – X-Ray window A-Position



Axis 350 – X-Ray window B-Position



Axis 350 – Hand Control

- Wired Hand Control (standard) with two connecting points on the operating table
- or
- Wireless Hand Control (optional) with RFID and Bluetooth technology

➤ Features:

- ✓ Colored backlit touch display 71 mm (2.8 inch)
- ✓ One Hand Control for multiple tables
- ✓ Individually programmable memory function
- ✓ Error messages and Log-file for service purposes
- ✓ Liquid protection: IPX7
- ✓ No need of extra charging station
- ✓ Motor drive control. No need for extra drive control unit



Item: **92523**

Wireless Hand control
(optional)

Axis 350 – Data Log on SD Card

Data Log supports failure analysis
Quick response on the failures



Axis 350 – Emergency Control

Manual emergency control system with foot pedal



Axis 350 – Foot Control

Possible Movements:

- ✓ Column up and down
- ✓ Trendelenburg / Reverse Trendelenburg
- ✓ Tilt right and left



Item: **92532**
Wired foot control
Axis 350



Item: **92533**
Wireless foot control
Axis 350

Axis 350 – Components

Axis 350 – High Version

Item: **92600H**
Axis 350
High Version




Axis 350 – Low Version

Item: **92600L**
Axis 350
Low Version




5th Wheel (optional)

Item: **92261**
Direction Control
Axis 350



Not motorized


Item: **92260**
Drive Unit
Axis 350



Motorized

Head Rest

Item: **92640**
Head rest




Intermediate Plates


Item: **92660**
Short seat plate



Item: **92631**
Long Intermediate plate




Item: **92630**
Short Intermediate plate




Leg Plates


Item: **92632**
Leg plates 2-part



Item: **92633**
Leg plate one-piece



Item: **92620**
Leg plate light



Axis 350 – Components

Connection Combinations




A - Position



B - Position

Combination	Connections						Length			
Combination 1	Head Plate 	+	Int. Plate Long 	+	Axis 350 	+	Leg Plates 	= 2180 mm		
Combination 2	Head Plate 	+	Int. Plate Short 	+	Axis 350 	+	Leg Plates 	= 2030 mm		
Combination 3	Leg Plates 		+	Axis 350 	+	Int. Plate Long 	+	Head Plate 	= 2180 mm	
Combination 4	Leg Plates 		+	Axis 350 	+	Int. Plate Short 	+	Head Plate 	= 2030 mm	
Combination 5	Head Plate 	+	Short Seat 	+	Axis 350 	+	Int. Plate Long 	+	Leg Plates 	= 2359 mm
Combination 6	Head Plate 	+	Short Seat 	+	Axis 350 	+	Int. Plate Short 	+	Leg Plates 	= 2209 mm

Basic Configuration



No	Item	Qty	Description
1	92600L 92600H	1	Axis 350 - Low Version (or) Axis 350 - High Version
2	92640	1	Head rest Axis 350
3	92630 92631	1	Short intermediate plate Axis 350 (or) Long intermediate plate Axis 350
4	92632 92633 92620	1	2-part Leg plates Axis 350 (or) Leg plate one-piece (or) Leg plate light
	91260	1	Packaging OT-Table WITHOUT ramp
	91261	1	Reusable exit ramp
Optional Accessories			
5	92660	1	Short seat plate
	92261	1	5th Wheel direction Control

Accessories for General Surgery



Item: **90321**

(or)



Item: **90323**

(or)



Item: **90635**

(or)



Item: **90636**

(+)



Item: **90506**

(or)



Item: **90310**

7



Item: **90520**

(or)



Item: **4646.05**

(+)



Item: **5.890-VE**

8



Item: **800-0057**

9



Item: **92021**


10



Item: **90500**

No	Item	Qty	Description
6	90321	1	Anaesthesia screen, rigid (+)
	90506	1	Radial setting clamp with fixing screw (or)
	90310		Rapid-action radial setting clamp
6	90323	1	Anaesthesia screen with central articulation (+)
	90506	1	Radial setting clamp with fixing screw (or)
	90310		Rapid-action radial setting clamp
6	90635	1	Anaesthesia screen, with extendable rods (+)
	90506	1	Radial setting clamp with fixing screw (or)
	90310		Rapid-action radial setting clamp
6	90636	1	Anaesthesia screen central art., ext. rods (+)
	90506	1	Radial setting clamp with fixing screw (or)
	90310		Rapid-action radial setting clamp
7	90520	2	Arm rest (or)
	4646.05		Arm support quick central fixation (+)
	5.890-VE		Fastening straps (1 pair)
8	800-0057	1	Restraint Strap
9	92021	1	Arm strap with fixing screw
10	90500	1	Infusion pole with clamp, height-adjustable
11	90332	1	Instrument tray with clamp, height-adjustable
12	90360	1	Accessories trolley (+)
	91300	1	Packaging for Accessories trolley
Dedicated Accessories			
13	92029	1	Universal Adapter for DORO® basic unit
14	90939	1	Table extension
15	90913	1	Hand surgery table with fixing screw
16	92006	1	X-ray cassette frame for Axis 350

Ophthalmology / ENT – Configuration A

1  Item: 92029	2  Item: 90456	3  Item: 90041	4  Item: 90003
--	--	---	--

	 <table> <tr> <th>No</th><th>Item</th><th>Qty</th><th>Description</th></tr> <tr> <td>1</td><td>92029</td><td>1</td><td>Adapter for DORO® basic unit, for Axis 350</td></tr> <tr> <td>2</td><td>90456</td><td>1</td><td>Hexagonal adapter with square adapter</td></tr> <tr> <td>3</td><td>90041</td><td>1</td><td>Head rest holder, with lateral adjustment</td></tr> <tr> <td>4</td><td>90003</td><td>1</td><td>Head rest (horseshoe) with pad</td></tr> </table>	No	Item	Qty	Description	1	92029	1	Adapter for DORO® basic unit, for Axis 350	2	90456	1	Hexagonal adapter with square adapter	3	90041	1	Head rest holder, with lateral adjustment	4	90003	1	Head rest (horseshoe) with pad
No	Item	Qty	Description																		
1	92029	1	Adapter for DORO® basic unit, for Axis 350																		
2	90456	1	Hexagonal adapter with square adapter																		
3	90041	1	Head rest holder, with lateral adjustment																		
4	90003	1	Head rest (horseshoe) with pad																		

Ophthalmology / ENT – Configuration B

1  Item: 92029	2  Item: 90045	3  Item: 1002.72A0	4  Item: 1004.75A0
--	--	---	--

No	Item	Qty	Description
1	92029	1	Adapter for DORO® basic unit, for Axis 350
2	90045	1	Horseshoe-headrest holder 3 joints
3	1002.72A0	1	Horseshoe-shaped head rest, two pieces
4	1004.75A0	1	Wrist hand support for ophtalmology

Abdominal – Configuration A

1 Item: 90520 Item: 4646.05 Item: 5.890-VE	2 Item: 90325 Item: 90506 Item: 90310	3 Item: 90330 Item: 90506 Item: 90310
---	--	--




No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90325 90506	1 1	Arm strap with fixing screw Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	90330 90506	2 2	Leg holder "Goepel" (1 piece) (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp

Abdominal – Configuration B

1 Item: 90520 Item: 4646.05 Item: 5.890-VE	2 Item: 90325 Item: 90506 Item: 90310	3 Adjustable Leg holder
---	--	---------------------------------------

No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90325 90506	1 1	Arm strap with fixing screw Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	800-0342-PR 800-0338-EU	1 2	Adjustable Leg holder for patients up to 182 kg (or) 800-0342-PL for patients up to 272 kg (or) 800-0364 for patients up to 454 kg (+) Springloc Clamp (order 2 pieces)



Bariatric, Laparoscopic – Configuration A

1  Item: 90520 (or)  Item: 4646.05 + 5.890-VE	2  Item: 90079 (+)  Item: 81501	3  Item: 800-0055 (+)  Item: 800-0338-EU	4  Item: 90939
---	---	---	--



No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90079 81501	1 2	Knee cradles (1 pair), with pads (+) Clamp, Knee cradles (order 2 pieces)
Optional			
3	800-0055 800-0338-EU	1 2	Nissen straps, 272 kg patient w. cap. (1 pair) (+) Springloc Clamp
4	90939	4	Table extension (1 piece),

Bariatric, Laparoscopic – Configuration B

1  Item: 90520  Item: 4646.05 + 5.890-VE	2  (+)  Adjustable Leg holder	3  Item: 800-0055 (+)  Item: 800-0338-EU	4  Item: 90939
---	--	---	--



No	Item	Qty	Description
1	90520	1	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	800-0342-PR 800-0338-EU	1 2	Adjustable Leg holder for patients up to 182 kg (or) 800-0342-PL for patients up to 272 kg (or) 800-0364 for patients up to 454 kg (+) Springloc Clamp (order 2 pieces)
3	800-0055 800-0338-EU	1 2	Nissen straps, 272 kg patient w. cap. (1 pair) (+) Springloc Clamp
4	90939	4	Table extension (1 piece),

Nephrectomy, Thorax

1



Item: 90520

(or)



Item: 4646.05
Item: 5.890-VE

(+)



2



Item: 90325

(+)



Item: 90506

(or)



Item: 90310

3



Item: 90351

(+)



Item: 90350

(+)

(or)



Item: 90539



Item: 90314

4



Item: 90352

(+)



Item: 90350

(+)

(or)



Item: 90539



Item: 90314

5



2

5



Item: 90609

(+)



Item: 90350

(+)


(or)



Item: 90539



Item: 90314



2

1

3



No	Item	Qty	Description
1	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90325 90506	1 1	Arm rest high, without clamp (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	90351 90350 90539	1 1 1	Side and back rest, with pad (+) Universal holder for side and back rest (+) Radial clamp with fixing screw, rotating (or) 90314 Rapid-action radial setting clamp
4	90352 90350 90539	1 1 1	Side and back supporting roll Universal holder for side and back rest (+) Radial clamp with fixing screw, rotating (or) 90314 Rapid-action radial setting clamp
5	90609 90350 90539	1 1 1	Sacrum/Pubis/Sternum support, with pad (+) Universal holder for side and back rest (+) Radial clamp with fixing screw, rotating (or) 90314 Rapid-action radial setting clamp

Rectal, Andrews, Genupectoral

1



Item: 90520

Item: 4646.05
Item: 5.890-VE

2



Item: 90277

Item: 90506

(or)

Item: 90310


3



Item: 90079


Item: 81501

4



Item: AP022

Item: QPHPD-ADULT 1



No	Item	Qty	Description
1	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
2	90277 90506	1 2	Abdominal rest, with pad, adjustable (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
3	90079 81501	1 2	Knee cradles (1 pair), with pads (+) Clamp, Knee cradles (order 2 pieces)
4	AP022	1	Gel Headrest for prone position (or) QPHPD-ADULT 1 – Head rest (single use)

Gynaecology / Urology – Configuration A

1

Item: 90330

Item: 90506
(or)
Item: 90310

2

Item: 9911202
(+)
Item: 9911012

3

Item: 90520
(or)
Item: 4646.05
Item: 5.890-VE



No	Item	Qty	Description
1	90330 90506	2 2	Leg holder "Goepel" (1 piece) (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
2	9911202 9911012	1 1	Tray TUR (+) Frame for Tray TUR
3	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)

Gynaecology / Urology – Configuration B

1

Adjustable Leg holder

2

Item: 9911202
(+)
Item: 9911012

3

Item: 90520
(or)
Item: 4646.05
Item: 5.890-VE


4

Item: 800-0074-S



No	Item	Qty	Description
1	800-0342-PR 800-0338-EU	1 2	Adjustable Leg holder for patients up to 182 kg (or) 800-0342-PL for patients up to 272 kg (or) 800-0364 for patients up to 454 kg (+) Springloc Clamp (order 2 pieces)
2	9911202 9911012	1 1	Tray TUR (+) Frame for Tray TUR
3	90520	2	Arm rest with clamp (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
Optional			
4	800-0074-S	1	Great White® Stirrup Dolly

Spine – Configuration


1  Item: NO 165	2  Item: NA 401	3  Item: 90520 Item: 4646.05 Item: 5.890-VE
---	---	--

4





Item: **92029** Item: **90456** Item: **90041** Item: **90003** Item: **QPHPD-ADULT 1** Item: **AP022**

5  Item: 800-0235	6  Item: 800-0074-SF
--	---



No	Item	Qty	Description
1	NO 165	1	Torso Pad "Onsernone"
2	NA 401	1	Leg support wave cushion 65x50x13cm
3	90520	2	Arm rest (or) 4646.05 Arm support quick central fixation (+) 5.890-VE Fastening straps (1 pair)
4	92029 90456 90041 90003	1 1 1 1	Adapter for DORO® basic unit, for Axis 350 (+) Hexagonal adapter with square adapter (+) Head rest holder, with lateral adjustment (+) Head rest (horseshoe) with pad (or) QPHPD-ADULT 1 - QView™ Head positioning device (single use) (or) AP022 Gel Headrest for prone position
Optional Accessories			
5	800-0235	1	SchureSpine Frame, Radiolucent
6	800-0074-SF	1	Schure Spine Frame Dolly

Neuro – Configuration A - Aluminium Skull Clamp System




No	Item	Qty	Description
1	92029	1	Adapter for DORO® basic unit, Axis 350 (or) 301110 Rail adapter (for tables not Schaerer)
2	3001-00	1	DORO® Base unit aluminum
3	3002-00	1	DORO® Swivel adaptor aluminum (or) DORO® Ball adapter, allows 360° rotation
4	1001.001	1	DORO® Skull clamp QR3 aluminum
5	3009-00	1	DORO® Horseshoe headrest for adults
6	3008-00	1	DORO® Horseshoe headrest for children
		1	Pins reusable (or) Pins disposable

Additional accessories for sitting position

7	3007-00 90506	1	DORO® Rail adapter for sitting position (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
---	------------------	---	---

Sitting Position



Pins reusable



Item: 3005-00 (Stainless Steel)
Item: 3005-50 (Titanium)
Reusable skull pins for **adults**
(3 pieces)



Item: 3004-00 (Stainless Steel)
Item: 3004-50 (Titanium)
Reusable skull pins for **children**
(3 pieces)

Pins disposable



Item: 1106.003
DORO® disposable skull pins, for **adults**,
(36 pieces)



Item: 3006-10
DORO® disposable skull pins, for **children**,
(36 pieces)

Pins disposable



Item: 3006-20
DORO® disposable skull pins, for **adults**,
(36 pieces)





Item: 3006-30
DORO® disposable skull pins, for **children**,
(36 pieces)

Neuro – Configuration B - Radiolucent Skull Clamp System


Item: 1101.020 (Complete Set)

DORO® Radiolucent Head Rest System compatible with all tables
(No adapter needed)



No	Item	Qty	Description
1	1101.021	1	DORO® LUCENT Basic support
2	1101.031	1	DORO® LUCENT Clamp intermediate piece
3	1101.026	1	DORO® LUCENT Intermediate piece
4	3033-00	1	DORO® Swivel adapter, radiolucent
5	3034-00	1	DORO® Skull clamp, radi. with Quick-Rail™
6		1	Pins reusable (or) Pins disposable
Additional accessories for sitting position			
7	3007-00 90506	1	DORO® Rail adapter for sitting position (+) Radial setting clamp with fixing screw (or) 90310 Rapid-action radial setting clamp
8	3001-00	1	DORO® Base unit aluminum

Sitting Position



Pins reusable



Item: 3005-00
Reusable skull pins
for adults (3 pieces)



Item: 3004-00
Reusable skull pins for
children (3 pieces)

Pins disposable



Item: 1106.003
DORO® disposable skull
pins, for adults,
(36 pieces)



Item: 3006-10
DORO® disposable
skull pins, for children,
(36 pieces)

Pins disposable



Item: 3006-20
DORO® disposable
skull pins, for adults,
(36 pieces)



Item: 3006-30
DORO® disposable
skull pins, for children,
(36 pieces)

ETS – Traction System

Applications

Femur




Tibia



Lateral Position



ETS - Femur

No	Images	Item	Description	No	Images	Item	Description
1		92690	Connection bar ETS	6		92693	Bar ETS (Changeable from left to right - order 2 pc.)
2		92691	Seat plate ETS	7		92694	Traction Unit ETS (Changeable from left to right - order 2 pc.)
3		90958	Offset holder for Perineal rod	8		92695	ETS Right Leg plate
4		90971	Perineal rod	9		92696	ETS Left Leg plate
5		90990	Perineal rod for Arthroscopy	10		90763	Extension shoe white (Changeable from left to right order 2 pc.)


ETS - Tibia Additional Items			
No	Images	Item	Description
11		92612	Traction Unit for Tibia ETS
12		92699	Foot plate
13		92698	Foot Plate Child
14		90193	Telescopic support

Optional Items			
No	Images	Item	Description
16		90520	Arm rest with rapid-action fixing clamp
17		90325 90506	Arm rest high (+) Radial setting clamp
18		90330 90506	Goepel (+) Radial setting clamp
19		81696	Side rail extension
20		90994	Extension bow holder

ETS - Lateral Position Additional Items			
No	Images	Item	Description
15		92611	ETS Counter Traction Unit for lateral position


ETS – Trolley			
No	Images	Item	Description
21		92613	ETS Trolley

ETS - 1 Leg Traction




No	Item	Qty	Description
1	92690	1	Connection bar ETS
2	92691	1	Seat plate ETS
3	90958	1	Offset holder for perineal rod (optional)
4	90971	1	Perineal rod
6	92693	1	Bar ETS
7	92694	1	Traction Unit ETS
8	92695	1	ETS Right Leg Plate
9	92696	1	ETS Left Leg Plate
10	90763	1	Extension shoe white
16	90520	2	Arm rest with rapid-action clamp
18	90330 90506	1	Goepel (+) Radial setting clamp
20	90994	1	Extension bow holder (optional)

ETS - 2 Leg Traction




No	Item	Qty	Description
1	92690	1	Connection bar ETS
2	92691	1	Seat plate ETS
3	90958	1	Offset holder for perineal rod (optional)
4	90971	1	Perineal rod
6	92693	2	Bar ETS
7	92694	2	Traction Unit ETS
8	92695	1	ETS Right Leg Plate
9	92696	1	ETS Left Leg Plate
10	90763	2	Extension shoe white (or)
16	90520	1	Arm rest with rapid-action clamp
20	90994	1	Extension bow holder (optional)

ETS – Tibia (in addition to Femur Traction)



No	Item	Qty	Description
11	92612	1	Traction Unit for Tibia ETS
12	92699	1	Foot plate
13	92698	1	Foot plate Child
14	90193	1	Telescopic support
18	90330 90506	1	Goepel (+) Radial setting clamp
Optional Accessories			
16	90520	2	Arm rest with rapid-action fixing clamp

ETS – Lateral Positioning (in addition to Femur Traction)



No	Item	Qty	Description
15	92611	1	ETS Counter Traction Unit for lateral position
Optional Accessories			
16	90520	1	Arm rest with rapid-action fixing clamp
17	90325 90506	1 1	Arm rest high (+) Radial setting clamp

Hand Surgery



1



2



(+)

Item: **90506**
(order 2 pc.)



1

No	Item	Qty	Description
1	90913	1	Hand surgery table, with fixing screw (or)
2	800-0040 90506	1 2	Hand surgery table, radiolucent Radial setting clamp

Humerus / Traction Device




1





No	Item	Qty	Description
1	9907038	1	Humerus and Forearm Positioning Device

Knee Arthroscopy



1

(or)

Item: **90357** Item: **90356**



(+)

Item: **90539**



(or)

Item: **90314**



2

(+)

Item: **90330**



(or)

Item: **90506**



(or)

Item: **90310**




1

2

No	Item	Qty	Description
1	90357	1	Leg holder arthroscopy (or)
	90356	1	Knee rest (+)
	90539	1	Radial clamp with fixing screw, rotating (or)
	90314	1	Rapid-action radial clamp, rotating
2	90330	1	Leg holder "Goepel" (1 piece) (+)
	90506	1	Radial setting clamp with fixing screw (or)
	90310	1	Rapid-action radial setting clamp



SSP - Schaerer Shoulder Plate - Configuration A


1  Item: 92625	2  Item: 90456	3  Item: 90041	4  Item: 90566 (or) Item: 90003
--	--	---	---



No	Item	Qty	Description
1	92625	1	Shoulder plate (SSP) Axis 350
2	90456	1	Adapter
3	90041	1	Head rest holder, with lateral adjustment
4	90566 90003	1	Helmet, radiolucent (or) Head rest (horseshoe) with pad

SSP - Schaerer Shoulder Plate - Configuration B

1  Item: 92625	2  Item: 92666	3  Item: 90899	4  Item: F-90023
--	--	---	--



No	Item	Qty	Description
1	92625	1	Shoulder plate (SSP) Axis 350
2	92666	1	Head rest holder
3	90899	1	Head rest
4	F-90023	1	Head pad (10 Pieces)

Axis 350 – One Table – Multiple Features



Data Log
on SD Card



Sliding 300mm



X-Ray Window
A-position: 1285mm



450 kg

Static Weight



Bluetooth wireless
Hand control unit



Emergency Panel



Pairing panel for
wireless
Hand control




Footpump



Quick click buttons



Double wheels 125mm
Cylinder Stamps



5th wheel
drive unit (*optional*)



Rubber Bellows

OP 830 - 850 **Basic Line**

Entry line Mobile Operating Table
suitable for many disciplines



Technical Specifications	OP 830	OP 850
Table load capacity	297 kg	350 kg
Table weight	275 kg	300 kg
Table length	2040 mm	2080 mm
Table width (without side rails)	520 mm	520 mm
Table range (without pads)	680 - 950 mm	650 - 1000 mm
Kidney-bridge elevator up	136 mm	128 mm
Head plate (up/down)	30° / 45°	30° / 45°
Back plate (up/down)	75° / 25°	80° / 45°
Leg plates (up/down)	5° / 90°	20° / 90°
Leg plates split	180°	180°
Trendelenburg / reverse	30° / 30°	30° / 30°
Lateral tilt (left/right)	20° / 20°	20° / 20°
Flex / reflex	205° / 105°	225° / 100°
Longitudinal shift electric	300 mm	300 mm
Wheels dimensions	120 mm	120 mm
Table top level reset	No	Yes
Exchange A&B position	No	Yes
Electro hydraulic function	Yes	Yes
Controls	Hand control	Hand contr. / Emerg. panel
Material (table top)	Powder coated iron	Stainless steel
Material (base cover)	Stainless steel	Stainless steel
Material (column cover)	Steel	Stainless steel
Pads height	60 mm	60 mm
Power supply	AC90V-264V / 50-60Hz	AC90V-264V / 50-60Hz

OP 830 / OP 850 - Ultrasound Welded Pads


OP 830 - Detachable plates (60 mm)

Memory foam

OP 850 - Detachable pads (60 mm)


Memory foam

Velcro fixation



OP 830 / OP 850 - Dimensions

Trendelenburg / Rev. Trendelenburg: **30°**




Lateral Tilt: **20°**




Back plate up: **80°**

Back plate down: **45°**




Head plate up: **30°**

Head plate down: **45°**




Leg plates up: **20°**

Leg plates down: **90°**




Leg plates abduction: **180°**




Leg plate abduction

OP 830 / OP 850 - X-Ray Dimensions

Position A




Position B



Segments

OP 830



OP 850



Range

OP 830



- Without pads: **680 – 950 mm**
- With pads: **740 - 1010 mm**
- Pads: **60 mm**

OP 850



- Without pads: **650 – 1000 mm**
- With pads: **710 - 1060 mm**
- Pads: **60 mm**

Safe Work Load Capacity

OP 830

297 kg



- Without pads: **680 – 950 mm**
- With pads: **740 - 1010 mm**
- Pads: **60 mm**

OP 850


350 kg



- Without pads: **650 – 1000 mm**
- With pads: **710 - 1060 mm**
- Pads: **60 mm**

Table Top Measurements

OP 830




OP 850



Sliding Top, Kidney Bridge

OP 830


300 mm
Electric sliding



Kidney bridge optional

OP 850

300 mm
Electric sliding



Kidney bridge optional

Controls

OP 830



Wired Hand control
(Wireless optional)



Emergency Panel on the column

OP 850



Wired Hand control
(Wireless optional)



Emergency Panel on the column

OP 830

Composition:

- Electro-hydraulic table
- Wired hand control (Wireless optional)
- Emergency control in the column
- Top of carbon steel
- Column and base stainless steel
- 60 mm memory standard pads
- Kidney bridge (optional)
- It can only be configured in position A
- Electric sliding 300mm
- Warranty 18 months

Included Components:

- Table (Seat and Back plates included)
- Head rest
- Leg plates
- Hand Control

Standard accessories included:

- 2 Arm rest (121000924)
- 1 Anaesthesia screen (121000247)
- 2 Knee crutch (121000295)
- 5 Round hole clamp (121000248)
- 2 Shoulder & Lateral support (121000296)



OP 850

Composition:

- Electro-hydraulic table
- Wired hand control (Wireless optional)
- Emergency control in the column
- Top of stainless steel
- Column and base stainless steel
- 60 mm memory standard pads
- Kidney bridge (optional)
- Can be configured in Position A and B
- Electric sliding 300mm
- "0" Level
- Warranty 18 months

Included Components:

- Table (Seat and Back plates included)
- Head rest
- Leg plates
- Hand Control

Standard accessories included:


- 2 Arm rest (121000924)
- 1 Anaesthesia screen (121000247)
- 2 Knee crutch (121000295)
- 5 Round hole clamp (121000248)
- 2 Shoulder & Lateral support (121000296)













OP 830 / OP 850 – The Different Models

Item	Description	Image
466901030	- OP 830 (<i>Head rest and Leg plates included</i>) - Standard accessories	
466901031	- OP 830 (<i>Head rest and Leg plates included</i>) - Standard accessories - Kidney bridge	
466901032	- OP 830 (<i>Head rest and Leg plates included</i>) - Standard accessories - Wireless Hand control	
466901033	- OP 830 (<i>Head rest and Leg plates included</i>) - Standard accessories - Kidney bridge - Wireless Hand control	
466901050	- OP 850 (<i>Head rest and Leg plates included</i>) - Standard accessories	
466901051	- OP 850 (<i>Head rest and Leg plates included</i>) - Standard accessories - Kidney bridge	
466901052	- OP 850 (<i>Head rest and Leg plates included</i>) - Standard accessories - Wireless Hand control	
466901053	- OP 850 (<i>Head rest and Leg plates included</i>) - Standard accessories - Kidney bridge - Wireless Hand control	

OP 830 / OP 850



1  Head plate	2  Clamp	3  Anaesthesia screen	4  Lateral support
5  Back plate	6  Arm rest	7  Kidney bridge (optional)	8  Body strap
9  Hand control	10  Knee crutch	11  Leg plates	12  Padding
13  Leg plate Abduction locking handle	14  Locking lever for Leg plate socket	15  Override control panel	16  AC power cable
17  Equipotential terminal	18  Seat plate	19  Table basement	20  Hand crank for Kidney bridge
21  Locking lever for Head plate socket	22  Connection for corded Hand control	23  Socket for AC mains	24  User manual Packing List



 Item: 121000350 Head rest + Item: 121000351 Adapter 	 Item: 121000347 Foot plate + Item: 121000347 Notched clamp (order 2 pieces)	 Item: 121000349 High position arm holder + Item: 121000248 Round hole clamp	 Item: 121000357 Drainage bowl	 Item: 121000355 Iliac support + Item: 121000248 Round hole clamp (order 2 pieces)
---	--	--	--	--

No	Item	Description
		Table
		OP 830 or OP 850
		<i>(The following accessories are included in the table)</i>
1	466901030 OP 830	2 Arm rest (121000924)
	466901050 OP 850	1 Anaesthesia screen (121000247)
		1 Body strap (121000251)
		2 Knee crutch (121000295)
		5 Round hole clamp (121000248)
		2 Shoulder & Lateral support (121000296)
1	825	1 Packaging
		Additional Acc. for Ophthalmology / ENT
2	121000350	Head rest
	121000351	Connection bracket for Neuro accessories
		Additional Acc. for Bariatric / Laparoscopic / Abdominal
3	121000347	Foot plate
	121000348	Notched clamp (order 2 pieces)
		Additional Acc. for Nephrectomy / Thorax / Kidney / Lateral
4	121000349	High position arm holder
	121000248	Round hole clamp
		Additional Acc. for Gynaecology / Urology
5	121000357	Drainage bowl
		Additional Acc. for Rectal / Andrews / Genupectoral
3	121000347	Foot plate
	121000348	Notched clamp (order 2 pieces)
6	121000355	Iliac support
	121000248	Round hole clamp



7

Item: **NO 165**
Torso Pad "Onsernone"

+

Item: **QPHPD-ADULT 1**
QViewTM Head positioning device



8

Item: **800-0235**
SchureSpine Frame, radiolucent

+

Item: **QPHPD-ADULT 1**
QViewTM Head positioning device



No	Item	Description
		Table
		OP 830 or OP 850 <i>(The following accessories are included in the table)</i>
1	466901030 OP 830	2 Arm rest (121000924) 1 Anaesthesia screen (121000247)
	466901050 OP 850	1 Body strap (121000251) 2 Knee crutch (121000295) 5 Round hole clamp (121000248) 2 Shoulder & Lateral support (121000296)
1	825	1 Packaging
		Additional Acc. for Spine Surgery (Option A)
7	NO 165 QPHPD-ADULT 1	Torso Pad "Onsernone" QViewTM Head positioning device (HPD)
		Additional Acc. for Spine Surgery (Option B)
8	800-0235 QPHPD-ADULT 1	SchureSpine Frame, radiolucent QViewTM Head positioning device (HPD)
		Additional Acc. for Shoulder Surgery
4	121000349 121000248	High position arm holder Round hole clamp
2	121000350 121000351	Head rest Connection bracket for Neuro accessories



Akrus Line

Mobile Surgical Tables and Chairs,
Breast Biopsy Chairs, Examination
Chairs for Ambulatory Surgery, Day
Surgery and Outpatient Surgery



Product Description

The mobile SC 5010 ES/HS surgical chair was developed in close cooperation with leading surgeons in the fields of Ophthalmology, ENT and Maxillo-facial surgery. The chair was designed to meet all the specific needs of surgeons working in the head, neck and shoulder area. It offers numerous ergonomic and functional user benefits.




- Operation and patient recovery on just one chair supports economic and efficient work flow
- Very comfortable and easy to patients as they do not have to move from one chair to another
- Chair height, back rest and seat cushion are individually electric
- Various head rests are available to offer unlimited patient head positioning
- Ample clear foot space underneath back rest
- Two types of column to choose from (ES model: 52 – 78 cm) or (HS model: 64 – 100 cm)

Dimensions and Weight


Length wheel base (backrest up)	760 mm
Width wheel base	580 mm
Width chair incl. Arm rests	700 mm
Width seat cushion	590 mm
Height chair back rest upright	1700 mm
Max. length chair, back rest horizontal	1750 mm
Wheel diameter	125 mm
Chair weight (depending on options)	114 kg

Column Range

AKRUS - SC 5010 ES



AKRUS - SC 5010 HS



Trendelenburg

AKRUS - SC 5010 ES




AKRUS - SC 5010 HS




Range Back Rest

AKRUS - SC 5010 ES




AKRUS - SC 5010 HS




Patient Weight Capacity

AKRUS - SC 5010 ES



AKRUS - SC 5010 HS





Features

- 1 – Headrest **electric adjustable**
- 2 – Height, Backrest, Seat and Leg Support: Individually **electric adjustable**
- 3 – Foot Control for castors on both sides
- 4 – Two height ranges possible: (52-78cm) or (64-100cm)
- 5 – Foot Switches
- 6 – Standard rails 25x10mm
- 7– Extending Foot Support

Applications

- Ophthalmology (ES)
- ENT Surgery (HS)
- Maxillofacial /Oral surgery (HS)
- Plastic esthetic surgery (HS)



Basic Configuration			
No	Item	Qty	Description
	242.000000 275.000000	1	Mobile surgical chair table SC 5010 ES (52 - 78 cm) (or) Mobile surgical chair table SC 5010 HS (54 - 100 cm)
1	241.030690 241.030660	1	Standard head rest (or) Head cushion, flat
5	242.012030 275.012020	1	Foot joystick, 5 segments (5 axes) (or) Foot switch, 2 segments (2 axes)
Optional Recommended Accessories			
8	4646.05	2	Arm rest lateral with 3 joints + (Item:5.890-VE Strap 1 pair)
9	240.030400	1	Arm support for surgeon on headrest, height adjustable

Product Description

The mobile SC 5010 SEK surgical chair table is developed in close cooperation with leading surgeons in various medical fields. The chair is designed to meet all the specific needs of surgeons. It offers numerous ergonomic and functional user benefits.

The infinitely adjustable backrest, the removable shoulder sections and various different adjustable head rests support any patient position.

- Removable shoulder sections
- Various head rests
- Versatile accessories
- Extremely low position
- Made in Germany



Applications

- Orthopaedics
- General Surgery
- Dermatology
- Plastic / Reconstructive Surgery

Dedicated Accessories

Item: 276.030648	Item: 276.030690	Item: 276.030700	Item: 276.030850

Dimensions and Weight

Length wheel base (backrest up)	760 mm
Width wheel base	580 mm
Width chair incl. Arm rests	700 mm
Width seat cushion	590 mm
Height chair back rest upright	1440 mm
Max. length chair, back rest horizontal	1960 mm
Wheel diameter	125 mm
Chair weight (depending on options)	95 kg

Features

- 1 – Height range: 68 – 98 cm
- 2 – Foot control for castors on both sides
- 3 – Standard Rails 25 x 10 mm
- 4 – Advanced battery concept
- 5 – Intuitive Hand control unit



Advanced Battery Management



Technical Data


Column Range



Trendelenburg



Range Back Rest




Patient Weight Capacity



Product Description

Patient chair for transporting and positioning patients for mammography examinations and stereotactic or tomosynthesis biopsies. The continuously adjustable backrest, the individually adjustable back segments and the height-adjustable headrest support any positioning for these treatments.

- Individually adjustable back segments / Arm rests for lateral positioning "decubitus position"
- Height-adjustable Head cushion
- In critical situations the shock position is available in seconds



Dedicated Accessories

Item: **277.032010**

Item: **277.030700**

Item: **9500300**


Applications

- Mammography
- Biopsy
- Tomosynthesis
- Stereotactic




Technical Data

Column Range



Trendelenburg



Range Back Rest



Patient Weight Capacity



Dimensions and Weight

Length wheel base (backrest up)	760 mm
Width wheel base	580 mm
Width chair incl. Arm rests	760 mm
Width seat cushion	620 mm
Height chair back rest upright	1440 mm
Max. length chair, back rest horizontal	1680 mm
Wheel diameter	125 mm
Chair weight (depending on options)	95 kg

Product Description



Technical Specifications

Base area W x L	620 x 680 mm
Weight chair	60 kg
Diameter dual castors	75 mm
Height adjustment	550 – 710 mm
Diameter seat cushion	420 x 420 mm
Thickness cushion	80 mm
Height adjustment arm	140 – 290 mm
Max. body weight	250 mm



The battery is located inside the foot switch housing. To recharge the battery connected the recharger to the jack in the back of the housing.

Accessories



Attachment Clamps

Item:90506

Radial setting clamp
with fixing screw,
rotating



- ➤ Radial adapter
- ⦿ ➤ Rotating
- ⚙ ➤ Gear Mechanism

Item:90310

Rapid-action radial setting clamp,
rotating



- ➤ Radial adapter
- ⦿ ➤ Rotating
- ⚙ ➤ Gear Mechanism

Item:90538

Fixing clamp with fixing screw,
non-rotating,
with square adapter



- ➤ Square adapter
- Non-Rotating

Item:90539

Radial setting clamp
with fixing screw, rotating,
with square adapter



- ➤ Square adapter
- ⦿ ➤ Rotating
- ⚙ ➤ Gear Mechanism

Item:90314

Rapid-action radial setting clamp,
rotating, with square adapter,



- ➤ Square adapter
- ➤ Radial adapter
- ⦿ ➤ Rotating
- ⚙ ➤ Gear Mechanism

Item:800-0228-EU

Simple Clamp



- Safety locks
- For accessories with flat blades

Item:800-0271-EU

Simple Clamp



- Accepts attachments with round poles up to 1.6 cm
- Attaches along the side rail without the need for notches

Item:800-0338-EU

Springloc Clamp for Leg holder
„Great White“




- Securely closes fittings
- Quick and easy release pull lever

Arm Rests

Item:90520

Arm support with rapid-action fixing clamp, 2 fixing strap with velcro fastener



- Viscoelastic pad
- Rapid-action fixing clamp
- Height adjustable
- Carrying capacity 15 kg
- 2 fixing straps with Velcro fastener
- Dim. Pad: 500 x 50 x 130 mm

Item:4646.05

Arm support, lateral

❖ Compatible strap:

- Item:5.890-VE (1 pair)



- Viscoelastic pad
- Anatomically shaped
- Articulated arm 460 mm
- Rails clamp
- Without fastening straps
- Dim. pad: 460 x 140 x 60 mm

Item:9906030

Arm support on ball joint with lever



- Can be swiveled via ball-and-socket joint by a lever system with one hand
- Height adjustable
- Completed with cushion dimension 510 x 140 mm
- Two fixation belts
- Frontal fast locking clamp.

Item:9906032

Arm support on ball joint



- Complete with an adjustment system with ball-and-socket joint to fix the support on all positions
- It is height adjustable with a run of 100 mm
- Complete with cushion 510 x 140 mm and thickness 40 mm
- Two fixation belts
- Frontal fast locking clamp

Item:9906019

Arm support for obese on ball joint



- Height adjustable for 150 mm
- Complete with frontal fast locking clamp
- Soft cushion
- Two fixation belts
- Dimensions: 500 x 248 mm.

Arm Rests

Item: **9906012**

Paediatric Arm support



- Completed with clamp for frontal and fast locking
- with a soft cushion
- Two fixation belts

Item: **90325**

Arm rest high

- ❖ Compatible clamps:
 - Art. 90506
 - Art. 90310



- For the top arm in lateral position
- With pad
- Height adjustable
- Carrying capacity 10 kg
- Fixing strap with Velcro fastener
- Dim. pad: 400 x 140 x 20 mm
- Also for positioning of the arm under the table

Item: **4670.03**

Arm support, vertical

- ❖ Compatible clamps:
 - Art. 90506
 - Art. 90310
- ❖ Compatible strap:
 - Item: **5.890-VE** (1 pair)



- Viscoelastic pad
- Anatomically shaped
- Articulated arm
- Without rails clamp
- Without fastening straps
- Dim. pad: 460 x 140 x 60 mm

Arm and Leg Guards

Toboggan Pads

Item: **508-0155** Toboggan Pad, Small
 Item: **508-0145** Toboggan Pad, Medium
 Item: **508-0144** Toboggan Pad, Large
 Item: **508-1399** Toboggan Pad, X-Large



- Available in four sizes: Small through X-Large
- Physician can lean against
- Simply slide under OT table pad to secure arms or legs (check pads fixation)
- For Axis 350 (not compatible with Axis and Arcus models)

Toboggan Arm/Leg Guards

- ❖ Available in four sizes:
 - Art. **800-0276** (27x61cm)
 - Art. **800-0082** (27x41cm)
 - Art. **800-0083** (20x41cm)
 - Art. **800-0084** (15x25cm)




- Simply slide under OR table pad to secure arms or legs
- Radiolucent
- For Axis 350 (not compatible with Axis and Arcus models)

Arm Holders

Item:NY 311

Arm holder "Grandson"



- In anodised aluminium
- Unbreakable
- Height 30 cm
- Width 23 cm
- With two universal rail holders

Item:NY 312

Arm holder "Grandson", short & adjustable version complete



- Short version
- Easy height adjustment directly at the operating table
- Anodised aluminium
- 2 x rail hangers with 2 fixing screws each
- Thermolabile soft padding with washable PU cover dark blue
- Height of top 204 mm
- Width plate 230 mm

Table Extensions

Item:92038

Table extension
for Arcus 501/601/701
(1 piece)



- For side rails 25 x 10 mm
- Dim. Pad: 470 x 175 x 210 mm

Item:90939

Table extension
for Axis 350/500/600/700
(1 piece)



- For side rails 25 x 10 mm
- Dim. Pad: 550 x 115 x 69 mm

Anaesthesia Screens

Item:90321

Anaesthesia screen, rigid,
with 2 drape clips.
without clamp.

- ❖ Compatible clamps:
 - Item:90506
 - Item:90310



- Height adjustable
- Tilttable
- Made of stainless steel
- Security ring

Item:90323

Anaesthesia screen with central
articulation

- ❖ Compatible clamps:
 - Item:90506
 - Item:90310



- Height adjustable
- Tilttable
- Central articulation
- Made of stainless steel

Item:90635

Anaesthesia screen,
rigid, with ext. rods
(2 pieces) 350 mm each.

- ❖ Compatible clamps:
 - Item:90506
 - Item:90310



- Height adjustable
- Tilttable
- Lateral extension bars
- Made of stainless steel
- Security ring

Item:90636

Anaesthesia screen with central
articulation, with ext. rods
(2 pieces) 350 mm each.

- ❖ Compatible clamps:
 - Item:90506
 - Item:90310



- Height adjustable
- Tilttable
- Central articulation
- Lateral extension bars
- Made of stainless steel

Anaesthesia screens

Item:9903003

Shaped anaesthesia screen

- ❖ Compatible clamp:
 - Art. 9901023



- Height adjustable
- Tilttable
- Rubber patient protection
- Security ring
- Complete with 2 lateral extension bar

Item:9903009

Flexible anaesthesia screen

- ❖ Compatible clamps:
 - (Order 2 pieces)
 - Art. 90506
 - Art. 90310
 - Art. 9901023



- Adjustable from 1400 to 1900mm
- Used with n. 2 clamps to be ordered separately

Item:3628.13

Anaesthesia tube holder



- Hose plate, folding arm 280 mm, column, base

Instrument tray

Item:90332

Instrument tray with fixing clamp



- Height adjustable
- With fixing clamp
- dimensions 520 x 330 mm

Thigh holder

Item:90326

Thigh holder



- With ball joint
- Without radial setting clamp
- With pad

Straps

Item:800-0057

Restraint Strap



- Amsco style
- All-rubber conductive strap
- Very durable
- Easy to clean and disinfect
- Mounts to surgical table with or without hooks
- Wide: 244 cm

Item:800-0107

Bariatric Restraint Strap



- Bariatric version
- Amsco style
- All-rubber conductive strap
- Very durable
- Easy to clean and disinfect
- Mounts to surgical table with or without hooks
- Wide: 353 cm

Item:9914012

Body/Leg strap with Velcro closure



- With clips for attaching onto the side rails
- With locking system, detachable in two sections
- Washable
- Width 100 mm.

Item:9914009

Body/Leg strap with clamps



- Width 100 mm
- Washable
- With clips for attaching onto the side rails adjustable in its length
- With anti fall system

Item:9914001

Body/Leg strap, width 60 mm

Item:9914008

Body/Leg strap for obese, width 150 mm



- Padded
- Washable belt
- Adjustable length
- With security and braking system
- Completed with clamps

Straps

Item:KM 700


Body strap "Lukmanier" 225 x 10cm



- Length 225 cm
- Width 10 cm
- Velcro length 55 cm
- Polyblac antistatic

Item:KG 710

Body strap "Klausen" 140 x 10cm



- Made of black polyester
- Washable
- Sterilisable at 134° C
- For leg attachment around leg plate for two legs

Item:800-0055

Nissen straps with 272 kg patient weight capacity (1 pair)


- ❖ Compatible clamps:
(order 2 pieces)
 - Art. 800-0338-EU



- Needed support during steep rev. Trendelenburg positioning
- Attaches to rail quickly and securely
- thick pad: 51 cm x 13 cm x 3 cm

Item:9914011

Single leg strap



- Washable
- With Velcro closure
- Suitable for divided leg section
- Qty advised 2

Item:KG 715

Leg plate strap "Umbrail"
168/60 x 10cm



- Made of black polyester
- 168 x 10 and 60 x 10 cm
- Washable
- Sterilisable at 134° C
- For fastening one leg around wide leg plate

Straps

Item: 9915001

Leather wristlet

Item: 9915003

Leather wristlet for paediatric patients

Item: 9915005


Leather wristlet for obese patients



- Washable
- With Velcro closure
- Padded
- Complete with clamp

Item: 9915002

Wrist support



- Made of washable cotton
- With Velcro closure
- Length: 100 mm
- Completed with fast-frontal locking clamp to be placed onto side rail to keep the arm parallel to the body

Item: 92021

Arm strap with fixing screw



- Dim: 590 x 105 mm
- With Velcro fastener
- Washable

Item: 9906801

Suspended arm strap



- Made of washable tissue
- 100 mm width
- Advised quantity 2 pieces
- To be attached to anaesthesia screen for the suspension of patient's arm

Item: Na 385

Arm suspension sleeve "Seelisberg" for head bow




- Used together with the head bow to secure the patient's arm in a hanging position
- The cuff is made of grey
- Padded imitation leather
- Has a window for infusion control

Lateral Supports

Item:90338

Side rest with fixing screw



- With pad
- Dim. pad: 220 x 100 x 30 mm


Item:90609

SPS sacrum- / pubis- / sternum support, with pad.




Item:90352

Side and back supporting roll with foam rubber pad



Item:90351

Side and back rest



Item:90350

Universal holder for side and back rest



Item:90538

Fixing clamp with fixing screw, **non-rotating**, with square adapter



Item:90539

Radial setting clamp with fixing screw, **rotating**, with square adapter




Item:90314

Rapid-action radial setting clamp, **rotating**, with square adapter,

Lateral Supports

Item:9908022
Universal Lateral Support


- ❖ Radial setting clamp included



- Made of stainless steel
- Swivelling, length and height adjustable with radial toothed locking system

❖ To be used with:

- Art. 9908001
- Art. 9908002
- Art. 9908004
- Art. 9908005
- Art. 9908006
- Art. 9908020
- Art. 9908026



Item:9908001

Sternal - Back Support

- Antistatic cushion
- Dim. Pad: 120 x 100 mm
- ❖ Requires: Art. 9908022



Item:9908002

Back Buttocks Support


- Antistatic cushion
- It can be mounted either horizontally or vertically
- Dim. pad: 280 x 140 mm
- ❖ Requires: Art. 9908022



Item:9908004

Chest Support


- Antistatic cushion
- Dim. pad: 170 x 200 mm
- ❖ Requires: Art. 9908022



Item:9908005

Sacral Support


- Antistatic cushion
- Can be mounted either horizontally or vertically
- Dim. Pad: 200 x 80 mm
- ❖ Requires: Art. 9908022



Item:9908006

Pubis Support

- Antistatic cushion
- Dim. Pad: 50 x 50 mm
- ❖ Requires: Art. 9908022



Item:9908020

Cylindrical Lateral Support

- Diameter: 80 mm
- Height: 300 mm
- ❖ Requires: Art. 9908022

Infusion Poles

Item:90500

Infusion pole with fixing clamp



- Extendable
- Can be tilted to maintain the vertical position
- Upper part with 4 bottle hooks
- Height-adjustment 850-1500 mm

Item:9916002

Infusion pole



- Complete with hooks and clamp
- Can be tilted to maintain the vertical position
- Upper part with 4 bottle hooks
- Height-adjustment 800-1500 mm

Side Rail Extension

Item:81696

Side rail extension



- Made of stainless steel
- For side rails 25 x 10 mm
- Length: 380 mm

Accessories Trolley

Item:90360

Accessories trolley
(without accessories)

Item:91300

Packing for accessories trolley




- With 4 side rails each 25x10 mm
- Equipped with 2 wire baskets
- Dimension: 580 x 500 mm
- Height: 1630 mm

X-Ray Cassette Telescopic Slide

Item: **92006**


X-ray cassette frame
for **Arcus 501/601/701**



- Telescopic slide: 750 – 1150 mm
- Dimensions:
 - Length: 340 mm
 - Width: 480 mm

Item: **90528**


X-ray cassette frame
for **Axis 500/600/700**



- Telescopic slide: 750 – 1150 mm
- Dimensions:
 - Length: 340 mm
 - Width: 430 mm

Item: **92670**

X-ray cassette frame
for **Axis 350**



- Telescopic slide: 750 – 1150 mm
- Dimensions:
 - Length: 460 mm
 - Width: 385 mm

X-Ray Cassette Holder

Item: **9922002** (1)


X-ray cassette holder for **all tables**
+

Item: **9922003** (2)

Cassette holder Frame



- Holding arm for positioning x-ray cassettes under radiotransparent table top
- Mounted on side rails by clamp (included)




- For inserting radiographic cassettes under the table
- Dimensions: 300 x 460 mm

Item: **90908**

X-ray cassette holder for **Axis 500**
with pad 69 mm

Item: **90909**

X-ray cassette holder for **Axis 600**
with pad 69 mm



- Consisting of: head-, back-, seat plate; 2 height equalisation pads for the leg plates
- The cassette channel is plugged-on instead of the standard pads
- Front side slide-in opening 450 mm

Head Rest Systems

Item:90003 (1)

Head rest (horseshoe)

+

Item:90041 (2)

Head rest holder, with lateral adjustment, with double articulated joint

+

Item:90456 (3)


Hexagonal adapter with square adapter

+

Item:92027 Adapter for **Arcus (4)**

Item:92028 Adapter for **Axis (4)**

Item:92029 Adapter for **Axis 350 (4)**



- For patients in supine position
- With antistatic cushion
- Adjustable
- ✓ Height
- ✓ Width
- ✓ Length

Item:9905039 (1)

Square Shaped Headrest

+

Item:90045 (2)

Headrest holder 3 joints

+

Item:92027 Adapter for **Arcus (3)**

Item:92028 Adapter for **Axis (3)**

Item:92029 Adapter for **Axis 350 (3)**



- For patients in supine position
- With antistatic cushion
- Adjustable
- ✓ Height
- ✓ Width
- ✓ Length

Item:9905011 (1)

Horse-Shoe Shaped Headrest

+

Item:90045 (2)

Headrest holder 3 joints

+

Item:92027 Adapter for **Arcus (3)**

Item:92028 Adapter for **Axis (3)**

Item:92029 Adapter for **Axis 350 (3)**



- For patients in supine and prone position
- With antistatic cushion
- Adjustable
- ✓ Height
- ✓ Width
- ✓ Length

Head Rest Systems

Item: **90001** (1)

Head rest

+

Item: **90041** (2)

Head rest holder, with lateral adjustment, with double articulated joint

+

Item: **90456** (3)


Hexagonal adapter with square adapter

+

Item: **92027** Adapter for **Arcus** (4)

Item: **92028** Adapter for **Axis** (4)

Item: **92029** Adapter for **Axis 350** (4)



- For patients in supine position
- With antistatic cushion
- Adjustable
- ✓ Height
- ✓ Width
- ✓ Length

Item: **1002.72A0** (1)

Horseshoe-shaped head rest, two pieces adjustable, with pad

+

Item: **90045** (2)

Horseshoe-headrest holder 3 joints

+

Item: **92027** Adapter for **Arcus** (3)

Item: **92028** Adapter for **Axis** (3)

Item: **92029** Adapter for **Axis 350** (3)



- For patients in supine position
- With antistatic cushion

- Adjustable
- ✓ Height
- ✓ Width
- ✓ Length

- Universal adapter for all Head Rest systems

Item: **1004.75A0** (4)

Wrist hand support for Ophthalmology

- ❖ Compatible Head rest:
Item: **1002.72A0**



- Arm support for the ophthalmologist

Adjustable Leg Holders

Item:800-0342-PR

Adjustable Leg holder for patients up to 182 kg

- ❖ Compatible clamps:
(order 2 pieces)
 - Item:800-0338-EU



- Rotates on three axes
- Lift assist allows easy control with minimal effort
- Flexible boot accommodates larger patients
- Patient weight capacity 182 kg
- Lithotomy range: -30° to 75°
- Adduction range: 9°
- Abduction range: +25°
- Externally rotates out: 20°

Item:800-0342-PL

Adjustable Leg holder for patients up to 272 kg

- ❖ Compatible clamps:
(order 2 pieces)
 - Item:800-0338-EU



- Rotates on three axes
- Lift assist allows easy control with minimal effort
- Flexible boot accommodates larger patients
- Patient weight capacity 272 kg
- Lithotomy range: -30° to 75°
- Adduction range: 9°
- Abduction range: +25°
- Externally rotates out: 20°

Item:800-0364

Great White HD 1000 Stirrups Boot System

- ❖ Compatible clamps:
(order 2 pieces)
 - Item:800-0338-EU



- Includes three pairs of Stirrup Boots: Premium, Platinum & Bariatric
- Patient weight capacity: 454 kg
- Lithotomy range: -25° to 75°
- Adduction range: 9°
- Abduction range: +25°
- Locks stirrup in place – stays locked until trigger is pulled
- Quick, easy connection
- Lithotomy indicators for accurate repositioning of legs

Leg Holders Dolly



Item:800-0074-S

Great White® Stirrup Dolly

- Made of stainless steel
- 4 swivel wheels
- Without Leg Holders

Leg Holders Wall Rack



Item:800-0284

Great White® Stirrup Wall Rack

- To be fixed to the wall for storing stirrups and clamps
- Without Leg Holders

Adjustable Leg Holders «Kids»

Item:800-0342-PUP

Pediatric E-Z Lift Stirrups for patients up to 73 kg

- ❖ Compatible clamps:
 - (order 2 pieces)
 - Item:800-0338-EU



- With exclusive axial rotation from the handle for superior patient positioning
- Boot adjusts with a single clamp which simplifies patient positioning
- Boot pads included
- Patient weight capacity 73 kg

Item:800-0342-KID

Pediatric E-Z Lift Stirrups for patients up to 113 kg

- ❖ Compatible clamps:
 - (order 2 pieces)
 - Item:800-0338-EU



- With exclusive axial rotation from the handle for superior patient positioning
- Boot adjusts with a single clamp which simplifies patient positioning
- Boot pads included
- Patient weight capacity 113 kg

Leg Holders «Goepel»

Item:90330

Leg holder «Goepel»
(1 piece)

- ❖ Compatible clamps:
 - Item:90506
 - Item:90310



- With pad
- without clamp
- With Velcro fixing straps (840 mm)
- Leg cradle adjustment through ball and socket joint
- Patient weight capacity 30 kg

Item:9902002

Knee crutch «Goepel»
(1 piece)

- ❖ Compatible clamp:
 - Art. 9901023



- With pad
- without clamp
- With Velcro fastener
- Leg cradle adjustment through ball and socket joint

Item:9902012

Paediatric Knee crutch «Goepel»
(1 piece)

- ❖ Compatible clamp:
 - Art. 9901023



- With pad
- with clamp
- With Velcro fastener
- Leg cradle adjustment through ball and socket joint
- Lengthwise tiltable, and laterally

Collecting Tray TUR

Item:90421

Collecting tray TUR
for **Axis 500/600**



- Made of nickel-chromium steel
- With strainer and drain hose
- With rapid action coupling
- Retractable
- Dim.: 440 x 310 x 180 mm.

Item:92037

Collecting tray TUR
for **Arcus 501/601**



- Made of nickel-chromium steel
- With strainer and drain hose
- With rapid action coupling
- Retractable
- Dim.: 345 x 265 x 140 mm.

Item:9911202 (1)

Tray TUR

+

Item:9911012 (2)

Frame for Tray TUR

❖ Compatibel with all tables



- Drainage SET
- With drain and recovery strainer
- Flexible tube to floor box

❖ Requires: Art. 9911012



- Fixed on side rails
- Made of stainless steel
- Adjustable
- Setting clamp included

❖ To be used with:
▪ Art. 9911202

Item:9911005

Floor Mobile Frame



- Made of stainless steel
- with 4 castors and handle
- Dimensions:
 - Height: 840 mm
 - Base Diameter: 385 mm

Abdominal Rest

Item:90277

Abdominal rest

- ❖ Compatible clamps:
 - (order 2 pieces)
 - Item:90506
 - Item:90310



- Adjustable
- With pad
- Without clamps (order 2 pieces)

Gynaecology handle

Item:9924021

Gynaecology handle
(order 2 pieces)



- To be fixed on the side guides
- Mounted fixing clamps

Knee Cradles

Item:90079

Knee cradles (1 pair)

- ❖ Compatible clamps:
 - (order 2 pieces)
 - Item:81501



- With pads
- Without clamps (order 2 pieces)
- Dim. Pads: 320 x 30 x 170 mm

Item:800-0285

Adjustable Foot Board/Table
Extension

- ❖ Compatible clamps:
 - (order 2 pieces)
 - Item:800-0338-EU



- Mounts vertically as a foot support when using reverse Trendelenburg
- Fits all OR tables and completely radiolucent
- Use as a table extension in the horizontal position to lengthen size of OR table
- Deluxe Pad included

Carbon-Fibre Plates 850 mm - Arcus

Item:90596

Carbon-fibre plate 850 mm
for **Arcus 501/701**

- ❖ On Seat interface
(instead of Seat plate)



- Unobstructed radiolucency
- Without side rails
- Dimensions: 850 x 520 x 130 mm
- Patient weight capacity: 220 kg

Item:90612

Carbon-fibre plate 850 mm
for **Arcus 501/601**

- ❖ On Leg interface
(instead of Leg plate)



- Unobstructed radiolucency
- Without side rails
- Dimensions: 850 x 520 x 130 mm
- Patient weight capacity: 90 kg

Carbon-Fibre Plates 850 mm - Axis

Item:90624

Carbon-fibre plate 850 mm
for **Axis 500/700**

- ❖ On Seat interface
(instead of Seat plate)



- Unobstructed radiolucency
- Without side rails
- Dimensions: 850 x 520 x 130 mm
- Patient weight capacity: 220 kg

Item:90625

Carbon-fibre plate 850 mm
for **Axis 500/600**

- ❖ On Leg interface
(instead of Leg plate)



- Unobstructed radiolucency
- Without side rails
- Dimensions: 850 x 520 x 130 mm
- Patient weight capacity: 90 kg

Carbon-Fibre Plate 1400 mm - Arcus

Item:90597

Carbon-fibre plate 1400 mm
for **Arcus 501/701**

❖ On Seat interface
(instead of Seat plate)



- Unobstructed radiolucency
- Without side rails
- Dim. : 1400 x 520 x 130 mm
- Patient weight capacity: 130 kg

Carbon-Fibre Plate 1400 mm - Axis

Item:90627

Carbon-fibre plate 1400 mm
for **Axis 500/700**

❖ On Seat interface
(instead of Seat plate)




- Unobstructed radiolucency
- Without side rails
- Dim. : 1400 x 520 x 130 mm
- Patient weight capacity: 130 kg

Child Surgery Plate

Item:90618

Child surgery plate
for **Arcus 501/701**
and **Axis 500/700**

❖ On Seat interface
(instead of Seat plate)



- With side rails
- Dim. Pad : 1060 x 400 x 69 mm
- Patient weight capacity: 50 kg

Additional Accessories to Carbon Plates

Item:90593

Side rail 25 x 10 mm
(1 piece)



- Slidable along the carbon plate
- Made in Inox 30 shot-blasted and electrolytically polished
- Dim. : 305 x 25 x 10 mm

Item:90536

Arm rest carbon-fibre for **Arcus**
(1 piece)

Item:90942

Arm rest carbon-fibre for **Axis**
(1 piece)




- Length of pad: 680 mm
- For insertion between pad and short (850 mm) or long (1400 mm) carbon-fibre plates
- Enables unobstructed X-raying from the heart to the fingers

Carbon Spine Frame (CSF)

Item:91600

Carbon Spine Frame (CSF)
for **Axis 500/700**
and **Arcus 501/701**

- Basic system includes:
 - ✓ Carbon frame
 - ✓ Head rest clamp including side rails short and mirror
 - ✓ 2 Cross bars
 - ✓ 2 Male chest pads
 - ✓ 2 Female chest pads
 - ✓ 2 Pelvic pads (unisex)
 - ✓ 3 Head rest QPHPD-ADULT 1
- The head rest clamp has the possibility to connect DORO® / MAYFIELD® skull clamps



Optional Accessories to Carbon Spine Frame

Item:91603

Side rail adapter **short**

Item:91604

Side rail adapter **long**



- In addition to standard side rails delivered with Carbon Spine Frame
- 2 Short side rails already included in the basic system

Item:4646.05

Arm support, lateral

❖ Compatible strap:

- Item:5.890-VE (1 pair)



- Viscoelastic pad
- Articulated arm 460 mm
- Including clamp
- Carrying capacity 10 kg
- Without fastening straps
- Dim. pad: 460 x 140 x 60 mm

Item:QPHPD-ADULT 1

QView TM Head Positioning Device (HPD) E.T. Tube Support and viewing mirror



- Single-use Head positioning device
- For prone position
- Mirror included

Item:AP022

Gel Prone Head Rest
28 X 23 X 14.5 cm

or

Item:AP023

Gel Prone Head Rest
23 X 19 X 13 cm



- Reusable Head positioning device
- For prone position
- Made of gel

Schure Spine Frame

Item:800-0235

SchureSpine Frame, Radiolucent



- Patient weight capacity: 227 kg
- Foam pads included
- Adjustable intraoperatively to the desired lordosis
- Tapered adjustable pads
- Radiolucent
- ✓ Schaerer Carbon Spine Frame compatible

Item:800-0074-SF

Schure Spine Frame Dolly
(without Spine Frame)



- For Schure Spine Frame

Torso Pads

Item:N0 165

Torso Pad "Onsernone"



- Horseshoe-shaped
- PU cover
- Washable

Item:N0 162

Torso Pad "Navegna"



- Box-shaped
- Dark blue PU cover
- Waterproof

Head Support

Item:9905035

Horse-Shoe Shaped Facial Support



- Width adjustable
- For patient in prone position
- Complete with:
- ✓ Cushions
- ✓ Supporting bar
- ✓ Clamp

DORO® Headrest Systems

Item:1001.100

DORO® QR3 Headrest System
(full set without pins)

❖ Requires:

Item:92027 Adapter for Arcus
Item:92028 Adapter for Axis
Item:92029 Adapter for Axis 350



- Full set without pins, consists of:
1. Item:1001.001 QR3 Skull Clamp
 2. Item:3002-00 Swivel Adaptor
 3. Item:3001-00 Adjustable Base Unit

Item:3003-100

DORO® QR3 Headrest System Teflon
(full set without pins)

❖ Requires:

Item:92027 Adapter for Arcus
Item:92028 Adapter for Axis
Item:92029 Adapter for Axis 350



- Full set without pins, consists of:
1. Item:1003-009 QR3 Skull Clamp
 2. Item:3002-009 Swivel Adaptor
 3. Item:3001-009 Adjustable Base Unit

Item:1101.020

DORO LUCENT® Headrest System
(full set without pins)



- Full set without pins, consists of:
1. Item:1101.021 Basic support
 2. Item:1101.031 Clamp interm. piece
 3. Item:1101.026 Intermediate piece
 4. Item:3033-00 Swivel adapter
 5. Item:3034-00 Skull clamp radiolucent with Quick-Rail™

Headrest Adapters

Item:92027

Adapter for DORO®, for Arcus

Item:92028

Adapter for DORO®, for Axis

Item:92029

Adapter for DORO®, for Axis 350



- For receiving:
- ✓ DORO® Head rest
 - ✓ MAYFIELD® Head rest
 - ✓ MIZUHO Head rest
 - ✓ Holders 900041, 90045, 90456

Item:3011-00

DORO® Side Rail OR Table Adaptor




- For receiving:
- ✓ DORO® Head rest
 - ✓ MAYFIELD® Head rest
 - ✓ MIZUHO Head rest

Item:3007-00

DORO® Crossbar Adaptor

- ❖ Compatible clamps
(order 2 pieces)
 - Item:90506



- To supplement DORO® QR3 Headrest System in a sitting position

Item:3009-00

DORO® Swivel Horseshoe Headrest
Adult w/ Extension Bar



- Universal non-invasive headrest designed for adults for prone or supine positioning
- Interchangeable gel pads

Item:3008-00

DORO® Swivel Horseshoe Headrest
Pediatric, with Extension Bar



- Universal non-invasive headrest designed for children for prone or supine positioning
- Interchangeable gel pads

Disposable Skull Pins			Reusable Skull Pins		
3006-00		DORO® Sterile Disposable Skull Pins, Stainless Steel, Adult (36 pcs)	3004-00		DORO® Reusable Skull Pins, Stainless Steel, Pediatric (3 pcs)
3006-10		DORO® Sterile Disposable Skull Pins, Stainless Steel, Pediatric (36 pcs)	3004-50		DORO® Reusable Skull Pins, Titanium, Pediatric (3 pcs)
3006-20		DORO® Sterile Disposable Skull Pins, Titanium, Adult (36 pcs)	3005-00		DORO® Reusable Skull Pins, Stainless Steel, Adult (3 pcs)
3006-30		DORO® Sterile Disposable Skull Pins, Titanium, Pediatric (36 pcs)	3005-50		DORO® Reusable Skull Pins, Titanium, Adult (3 pcs)
3006-50		DORO® Sterile Disposable Skull Pins, Stainless Steel, Adult (36 pcs)			
1106.003		DORO LUCENT® Sterile Disposable Skull Pins, Adult (3 pcs)			
1106.103		DORO LUCENT® Sterile Disposable Skull Pins, Adult (9 pcs)			

MTS Modular Traction System

- Available for:
- ✓ Axis 500/700
- ✓ Arcus 501/701

One surgical table – Multiple disciplines!

- All in one system for all surgical disciplines without change of seat plate
- Radiolucent carbon seat plate

Applications

General Surgery



Femur



Tibia



Hip



Lateral Position



Shoulder



Interface for
standard accessories

Interface for
perineal rod

Carbon-fibre
seat plate

Interface for
traction bar


Interface for
supplement bar



- MTS Trolley
- Easy storage
- Measurements: 625 x 890 x 1610 mm



MTS - Accessories – Arcus 501/701			MTS - Accessories – Axis 500/700		
No	Item	Description	No	Item	Description
1	90658	Seat Plate ADVANCED, Arcus 501/701	1	90962	Seat Plate ADVANCED, Axis 500/700
2	90661	Supplement bar, left, Arcus 501/701	2	90989	Supplement bar, left, Axis 500/700
3	90662	Supplement bar, right, Arcus 501/701	3	90988	Supplement bar, right, Axis 500/700
4	90971	Perineal rod	4	90971	Perineal rod
5	90990	Perineal rod for arthroscopy	5	90990	Perineal rod for arthroscopy
6	90958	Offset holder	6	90958	Offset holder
7	90983	Additional joint	7	90983	Additional joint
8	90981	Traction bar Metal	8	90981	Traction bar Metal
9	90974	Traction bar ADVANCED	9	90974	Traction bar ADVANCED
10	90978	Multiclamp	10	90978	Multiclamp
11	90928	Multiclamp slide	11	90928	Multiclamp slide
12	90980	Traction unit	12	90980	Traction unit
13	90895	Traction unit, vertical shift	13	90895	Traction unit, vertical shift
14	90763	Extension shoe MTS, white (1 piece)	14	90763	Extension shoe MTS, white (1 piece)
15	RO.0162.2018	Traction shoe with side rail interface	15	RO.0162.2018	Traction shoe with side rail interface
16	RO.0148.2018	Pad set incl. fixation straps (10 pieces)	16	RO.0148.2018	Pad set incl. fixation straps (10 pieces)
18	90994	Extension bow holder	18	90994	Extension bow holder
19	90780	Central rail adapter	19	90780	Central rail adapter
20	90897	Radial setting clamp with star screw	20	90897	Radial setting clamp with star screw
21	90759	Holder for knee rest	21	90759	Holder for knee rest
22	90896	Adapter for lower leg surgery	22	90896	Adapter for lower leg surgery
23	92699	Foot plate Adult (92698 for child)	23	92699	Foot plate Adult (92698 for child)
24	90193	Telescopic support	24	90193	Telescopic support
25	90688	Head rest system, Arcus 501/701	25	90762	Head rest system, Axis 500/700
26	F-90023	Head Pad for MTS (10 pieces)	26	F-90023	Head Pad for MTS (10 pieces)
27	90745	Extension plate, Arcus 501/701	27	90948	Extension plate, Axis 500/700
28	90677	Lateral ext. plate, left, Arcus 501/701	28	90925	Lateral ext. plate, left, Axis 500/700
29	90675	Lateral ext. plate, right, Arcus 501/701	29	90924	Lateral ext. plate, right, Axis 500/700
30	92611	Counter Traction Unit for lateral pos.	30	92611	Counter Traction Unit for lateral pos.
31	90670	Gyn. pad set, Arcus 501/701	31	90931	Gyn. pad set, Axis 500/700
32	90673	Gyn. cut-out plate, Arcus 501/701	32	90959	Gyn. cut-out plate, Axis 500/700
33	90957	MTS Trolley	33	90957	MTS Trolley
34	90738	Bar MIS (1 piece)	34	90738	Bar MIS (1 piece)
35	90775	MIS Trolley	35	90775	MIS Trolley



ETS – Accessories - Axis 350		
Nr	Item	Description
1	92690	ETS Connection Bar
2	92691	ETS Seat Plate
3	90958	Offset holder, for perineal rod
4	90971	Perineal rod
5	90990	Perineal rod for arthroscopy
6	92693	ETS Bar (1 piece)
7	92694	ETS Traction Unit
8	92695	ETS Leg Plate right
9	92696	ETS Leg Plate left
10	90763	Extension shoe, white (1 piece).
11	90994	Extension bow holder, for Kirschner-wire or pin
12	92612	ETS Traction Unit for Tibia
13	92699	Foot plate adult, length: 300 mm, width: 100 mm
14	92698	Foot plate child, length: 240 mm, width: 80 mm
15	90193	Telescopic support
16	92611	ETS Counter Traction Unit for lateral position
17	92613	ETS Trolley

Hand Surgery Tables

Item:90260

Hand surgery table for **Arcus**

Item:90913

Hand surgery table for **Axis**



- With fixing screw
- Without floor support
- With pad
- Dim. Pad Arcus:
 - 715 x 330 x 44 mm
- Dim. Pad Axis:
 - 710 x 330 x 69 mm

Item:800-0027

Carbon Fiber Minor Procedure Table

- ❖ Compatible clamps:
 - (order 2 pieces)
 - Item:800-0248-EU



- Adjusts up or down
- Table measures (38cm x 76cm)
- 2 required clamps separately
- Without floor support
- Deluxe foam pad included

Item:800-0040

Hourglass procedure table

- ❖ Compatible clamps:
 - (order 2 pieces)
 - Item:90506



- Radiolucent
- With pad
- Dimension: 380 x 840 x 580 mm

Shoulder Rests

Item:90334

Shoulder rests (1 pair)

- ❖ Compatible clamps:
 - (order 2 pieces)
 - Art. 90538



- Height and width adjustable
- With antistatic pads
- Without clamps
- Dim. pad: 220 x 100 x 30mm

Knee Arthroscopy Supports

Item:90356

Knee rest

- ❖ Compatible clamps:
 - Item:90539
 - Item:90314



- Without radial setting clamp
- With foam rubber pad
- Pad diameter 80 mm
- Pad length: 220 mm

Item:90357

Leg holder arthroscopy

- ❖ Compatible clamps:
 - Item:90539
 - Item:90314



- Without radial setting clamp
- With 2 adjustable thigh supports
- Adjustable height
- Radial swivelling
- With pad

Item:800-0065

Schure Grip Leg Holder

- ❖ Compatible clamps:
 - Item:90506
 - Item:90310



- Dual independent braces freely rotate 360° and spread horizontally from 8-28 cm to fit any size leg at any angle
- Use also for ACL procedures
- Made of stainless steel
- Leg support: 31 cm x 4 cm
- Braces: 20 cm
- Mounting post: 41 cm x 1.6 cm
- Schure Socket XL sold separately

Item:9907075

Foldable knee arthroscopy device



- Reclinable
- Spreadable
- Tiltable lateral supports
- Adjustable lateral supports
- Complete with padding and clamp

Item:276.032700

Knee clamp, electric



- The battery is located inside the foot switch housing
- To recharge the battery connected the recharger to the jack in the back of the housing



Knee Positioner

Item:800-0141

Total knee replacement positioner

- ❖ Compatible clamp:
 - Item:800-0271-EU



- Locks in any position with a simple turn of the wrist
- Smooth sliding motion for effortless leg positioning
- Well leg lies flat – eliminating the need for well-leg positioners, sandbags and gel pads
- Easy-to-use and completely rigid once locked
- All components are autoclavable
- Ball joint allows the surgeon to position the operative leg in extension, flexion, tilt and rotation

Hip Prothesis Positioning Device

Item:9908027

Hip Prothesis Positioning Device



- It is made of 2 height adjustable iliac crest tampon supports and 1 orientable and height adjustable sacral support
- Complete with quick-lock clamps to be installed on the side rails

Humerus / Traction Device

Item:9907038

Humerus and Forearm Positioning Device

- ❖ All items included



Item:9909013

Height adjustable plate for humerus surgery



- Radiolucent
- Orientable on ball-and-socket joint
- To use with patient in prone position
- Complete with cushion and clamp for side rails
- Height adjust. range: 150mm

Shoulder Plates for Arcus, Axis and Axis 350

Item:92625


Shoulder plate (SSP) Axis 350

Item:92645

Shoulder plate (SSP) Axis 500-700

Item:92655

Shoulder plate (SSP) Arcus 501-701



Configuration A

Item:92666 (1)

Head rest holder

+


Item:90899 (2)

Head rest

+

Item:F-90023 (3)

Head pad
(10 pieces)



Configuration B

Item:90456 (1)

Hexagonal adapter

+

Item:90041 (2)

Head rest holder, with lateral adjustment, with double articulated joint

+


Item:90003 (3)

Head rest (Horseshoe)

or

Item:90566 (4)

Helmet, radiolucent., with pad and fixing strap



Shoulder Plate for all Surgical Tables

Item: **800-0142-EU**

E-Z Lift beach chair (EU), patient weight capacity: 227 kg. Including Arm rest.



- Double-ball joint head rest
- Adjustable to any OR table
- Patient weight capacity: 227 kg
- Articulation range: 0°–90°

Shoulder Positioning Accessories

Item: **800-0169**

Lateral shoulder holder

- ❖ Compatible clamp:
 - Item: **800-0271-EU**



- Weightless
- Internal and external rotation can be easily controlled by the surgeon

Item: **800-0275**

Elbow surgical positioner

- ❖ Compatible clamp:
 - Item: **800-0271-EU**




- Intraoperative positioning for dorsal elbow surgery
- Quick release allows for intra-operative range of motion verification

Item: **800-0280-EU**

Schüre Loc™ XPS

- ❖ Additional optional accessories:
 - Item: **790-0039**
 - Item: **790-0040**
 - Item: **790-0056**
 - Item: **790-0057**
 - Item: **790-0077**



- Provides rigid fixation
- Unlimited intraoperative range
- Rotator cuff
- Replacement surgical procedures
- Trimano like movements

Schure Loc™ XPS Portable arm and hand positioner

- Trimano like movements and advanced features
- Provides rigid fixation
- Unlimited intraoperative range
- Rotator cuff
- Replacement surgical procedures

Item:800-0280-EU
Schure Loc™ XPS




Item:790-0281
Disp. Schure Loc™ XPS Sterile Kit,
10/cs
Includes:
✓ Sterile Forearm Wrap
✓ Sterile Self-Adherent Coban™
✓ Sterile Drape (fits whole arm)



Item:790-0039
Schure Loc™
XPS Adapter




Item:790-0056
Schure Loc™ XPS Armboard



Item:790-0040
Schure Loc™ XPS Forearm Attachment



Item:790-0057
Schure Loc™ XPS Counter Traction Post



Item:790-0075
Schure Loc™ XPS Elbow/Ankle Positioner



Item:790-0077
Schure Loc™ XPS TKR Positioner



Head Cushions

Item:NA 200

Head cushion with dent
20 x 20 x 6 cm



- It optimally supports the patient's head in the supine position
- The cushion is made of imitation leather

Item:NA 201

Cushion with dent and fastening compartments
20 x 20 x 6 cm



- Supports the patient's head in the supine position
- Made of imitation leather
- Equipped with fastening pockets on both sides

Item:NA 202

Cushion with dent and fastening compartments
30 x 20 x 6 cm



- Supports the patient's head in the supine position
- Made of imitation leather
- Equipped with fastening pockets on both sides

Item:NA 220

Cushion with dent, wedge-shaped
25 x 25 x 4-8cm



- It optimally supports the patient's head in the supine position
- Made of imitation leather

Shoulder Cushions

Item:NO 180

Shoulder support cushion "Muggio"



- For Trendelenburg positioning
- PU cover
- Washable
- With side compartments for tube holder

Item:NO 180-1

Shoulder support cushion "Muggio" with support plate



- For Trendelenburg positioning
- With support plate
- PU cover
- Washable
- Anti-slip coating
- With side compartments for tube holder

Chest Cushion

Item: NO 170

Chest support pad «Maggia» 1-piece



- For positioning the upper body in the prone position
- Made of viscoelastic SAF
- Has a PU cover in dark blue with welded seams
- Waterproof and washable

Leg Cushions

Item: NO 425

Three-quarter roll cushion «Generoso»
50 x 20 x 12 cm



- For two legs
- Dark blue PU cover
- Thermolabile
- Dimensions: 50 x 20 x 12 cm

Item: NO 426

Three-quarter roll cushion «Agnò»
25 x 20 x 12 cm



- For one leg
- Dark blue PU cover
- Thermolabile
- Dimensions: 50 x 20 x 12 cm

Item: NA 401

Leg support wave cushion
for **both legs**



- Suitable for supporting both legs
- Covered with imitation leather
- Dimensions: 65 x 50 x 13 cm

Item: NA 402

Leg support wave cushion
for **one leg**



- Suitable for supporting one leg
- Covered with imitation leather
- Dimensions: 65 x 22 x 13 cm

Tunnel Cushions

Item: NO 460

Tunnel cushion "Leventina"
2-Parts



- Thermolabile SAF
- With dark blue PU cover
- Waterproof
- Length 70 cm.

Item: NO 465

Module tunnel cushion "Airolo"



- 4-piece
- Dark blue PU cover
- Velcro fasteners
- Total length 80 cm.

Item: NA 460

Tunnel cushion



- Suitable for supporting legs in lateral position
- Covered with imitation leather
- Dimensions: 80 x 35 x 18-20 cm
- Radiolucent

Heel Cushions

Item: NO 455

Heel cushion pair





- Pair of heel pads thermoactive
- With dark blue PU cover
- Sewn together
- Foldable
- Washable
- Dimensions: 270 x 200 x 200 mm

Item: NO 456


Heel support splint cushion "Vedeggio"
(per piece)








- Sold by the piece
- PU cover in dark blue
- Waterproof
- Washable
- Dimensions: 43 x 15cm





Description	Image	Item	Dimensions (cm)
Operating Table Pad		AP301	20 X 50 X 1
		AP302	99 X 50 X 1
		AP303	50 X 50 X 1
		AP304	50 X 25 X 1
		AP305	10 X 38 X 1
		AP306	42 X 15 X 1
		AP307	10 X 61 X 1
		AP308	50 X 45 X 1
		AP401	193 X 86 X 1
		AP412-1	46 X 33 X 1
		AP413	22 X 17 X 1
		AP416	75 X 25 X 2
		AP417	10 X 15 X 2
		AP418	16 X 8 X 0.6
		AP423	70 X 50 X 1.6
		AP555	66 X 15 X 1.3
		AP701	97 X 64 X 1
Operating Table Pad		AP401	193 X 86 X 1
		AP402	183 X 50 X 1
		AP403	117 X 50 X 1
		AP404	160 X 40 X 1
		AP053	31 X 66 X 1.27
Operating Table Pad		AP410W/C	100 X 50 X 1
		AP411W/C	50 X 50 X 1
		AP409W/C	117 X 50 X 1
3 pieces Operating Table Pad		AP304	50 X 25 X 1
		AP410W/C	100 X 50 X 1
		AP411W/C	50 X 50 X 1
Embossed Operating Table Pad		AP061	50 X 50 X 1.6
		AP062	115 X 50 X 1.6
		AP063	40 X 23 X 1.6
		AP064	91 X 45 X 1.6



Description	Image	Item	Dimensions (cm)
Closed Head Ring		AP001	9 X 4 X 2
		AP002	14 X 6 X 3.5
		AP003	20 X 7.5 X 4.5
		AP004	20 X 7.5 X 7.5
Open Head Ring		AP011	9 X 4 X 2
		AP012	14 X 6 X 3.5
		AP013	20 X 7.5 X 4.5
		AP014	20 X 7.5 X 7.5
Open Head Ring		AP011/C	9 X 4 X 2
		AP012/C	14 X 6 X 3.5
		AP013/C	20 X 7.5 X 4.5
Prone Head Rest		AP022	28 X 23 X 14.5
		AP023	23 X 19 X 13
Head Pad		AP110	27 X 23 X 7
Supine Head Rest		AP112	27 X 24 X 9
Contoured Head Rest		AP081	19 X 19 X 6
Face Pad		AP021	16 X 6 X 1.6




Description	Image	Item	Dimensions (cm)
Hand / Wrist Supporter		AP005	28 X 8.3 X 6
Arm Board Pad		AP071	40 X 12.5 X 2
		AP072	60 X 12.5 X 2
Arm Protector		AP311	46 X 21 X 0.8
Arm Protector		AP310	38 X 12 X 0.8
Arm Board Pad with Foam		AP073	40 X 12.5 X 6
		AP074	60 X 12.5 X 6
		AP075	50 X 12.5 X 6
Contoured Arm Board Pad (Pair)		AP077	52 X 15 X 5.5
Heel Support Pad (Pair)		AP031	20 X 11 X 7
Heel Support Pad (Pair)		AP032	17 X 10.5 X 5

Description	Image	Item	Dimensions (cm)
Comfort Heel Cups with Strap (Pair)		AP033	17 X 18 X 1.3
Comfort Heel Cups with Strap (Pair)		AP034	17 X 23 X 1.3
Crutch Stirrup Pad (Pair)		AP801	31 X 23 X 1
		AP802	40 X 25 X 1
		AP803	52 X 30 X 1
		AP804	55 X 52 X 1
Boot Stirrup Pad		AP901	23 X 28 X 1.3
Traction Boot Pad Cover		AP902	37 X 24 X 0.5
Lithotomy Pole Covers (Pair)		AP903	71 X 12 X 0.8
Ankle Strap Stirrups (Pair)		AP904	31 X 4 X 1
		AP905	44 X 45 X 1
Gel Fastening Strap (Pair)		FB-01	33 X 7 X 1.5
		FB-02	36 X 7 X 1.5

Description	Image	Item	Dimensions (cm)
Chest Rolls (Adult)		AP041	30 X 15 X 10
		AP042	40 X 15 X 10
		AP043	50 X 15 X 10
		AP047	50 X 18 X 15
Chest Rolls with flat Bottom		AP047-1	15 X 5 X 5
		AP047-2	20 X 5 X 5
		AP047-3	25 X 5 X 5
		AP047-4	30 X 5 X 5
Chest Rolls (Paediatric)		AP048-1	15 X 7.5 X 7.5
		AP048-2	20 X 7.5 X 7.5
		AP048-3	25 X 7.5 X 7.5
		AP048-4	30 X 7.5 X 7.5
Contoured Chest Rolls		AP049-1	30 X 16 X 13/9
		AP049-2	40 X 16 X 13/9
		AP049-3	50 X 16 X 13/9
Chest Positioning Set		AP050-1	71/58 X 20 X 7.6
		AP050-2	71/58 X 20 X 7.6
Dome Positioner		AP051	35 X 15 X 7.5
Universal Patient Positioner		AP091	16 X 8 X 4.5
		AP093	33 X 12 X 6
		AP094	36 X 12.5 X 7
		AP095	25 X 12.5 X 5
		AP096	38 X 12.5 X 6.3
		AP097	50 X 12.5 X 7

Description	Image	Item	Dimensions (cm)
Graduated Thoracic Positioner		AP052	23 X 36 X 3.8
Bassinet Pad		AP053	30.5 X 66 X 1.27
Fracture Table Pad		AP054	14 X 20 X 1.3
Fracture Table Pad		AP055	18 X 18 X 1.3
Fracture Table Pad		AP056	33 X 30 X 0.8
Operating Table Pad with Foam		AP419	70 X 20 X 8
		AP420	30 X 15 X 10
		AP421	50 X 15 X 10
		AP422	60 X 15 X 10
FractureTable Post Pad		AP556	15 X 13 X 2.5
Patient Positioning Belt		AP601	50.5 X 8.5 X 0.8
Slanting Shaped Pad		AP102-3	30 X 22 X 5/1
		AP102-6	50 X 12 X 10/4
Sacral Protector		AP103-1	40 X 38 X 2.5
Thyroid Pad (does not include the head ring)		AP104-1	30 X 40 X 12
Thyroid Pad (does not include the head ring)		AP104-2	48 X 39 X 9

Description	Image	Item	Dimensions (cm)
Lateral Positioner		AP101	79 X 50 X 10
Tunnel Pad		AP103	41 X 65 X 2 (PAD) 41 X 65 X 21.5 (FRAME)
Angio Catheter Table		AP105	50 X 100 X 25
Relton Frame Pad Set (4 pieces)		AP111	19 X 17.5 X 10.5
Wilson Frame Pad		AP210	75 X 25.5 X 2.5
		AP211	75 X 25.5 X 1
Wilson Frame Pad Set		AP212	72.5 x 17.5 x 2
		AP213	60 X 49 X 15.5
Wilson Frame Set		AP214	67 X 16 X 2 (2pcs) 13 X 15.5 X 2 (2pcs) 15.5 X 15 X 2 (2pcs)
		AP215	47 X 67 X 2



Head office Switzerland




Branch office USA



- ✓ **Design innovation**
Copied often – never reached
- ✓ **Material innovation**
Carbon accessories
- ✓ **Engineering innovation**
Best-in-class features and versatility
- ✓ **Ergonomic innovation**
Tables and accessories designed for ease of use
- ✓ **Safety and hygiene innovation**
Patient safety and hygiene features

Partner network in over 75 countries





Schaerer Medical AG
Erlenauweg 17
CH-3110 Münsingen
Tel. +41 31 720 22 00
info@schaerermedical.ch
www.schaerermedical.com