



MCB Series

MCB - 14W MCB14W Single Crank Adjusment

Manual Bed

Features

MCB14W Single Crank Adjusment Manual Bed

Side-rail-free beds with single crank adjustment of backrest, fixed height. It's suitable to meet basic care needs of patients in need. The bed to be controlled by caregivers for short and long term.

Suggested area of use : Home Care, Hospitals, Nursing

Patients can experince an ideal sitting position with adjustment of backrest to 0-70°, therefore activities like eating, using bedpan, watching tv, reading boo or using laptop to be done comfortably.

Side Rails: Developed as side-rail-free, for those who do not have risk of falling and have auto-control capability providing cost effective solutions. Main target users are patients who have demantia, alzheimer, panic attack etc. Please consult your doctor or an expert before order the bed.

Head and Foot Panels: They are designed as detachable for emergency intervenes as well as making up the bed. Wood/Metal type panels developed for being economic and practical. It easily fit on sockets placed at four corner of the bed. Metal part is stainless steel, and wood part is MDF covered wood. Height of the panels is 33,5 cm from mattress platform. Metal parts of the bed made of smooth finishing for injury risks.

Mattress Platform: 2 sectional beds, one is adjustable backrest to 0-70°, another one is fixed section. Mattress platform has wide lying area as 85cm width and 190cm lenght for a comfortable recovery.

Painting and Coating: Construction and frame surfaces are coated with epoxy powder paint at high degree ovens. Durability, resistance to scratching and smooth surface at high quality is under guarantee. Therefore, It keeps the look like the first day. Eco-friendly smooth surfaces help you take care of your beds at ease.

Construction (Body) : Made of steel construction. With gas metal arch welding technique, It provides strong and durable body.

Additional Option: Legs and boards of our beds are demountable, therefore It ensures maximum use of space during storage and transportation. There is space under beds needed for placing lifts and overbed tables.

Adjustment Control : Adjustments made via hand crank control. Beds equiped with mechanism that helps caregivers adjusting the bed smoothly and in silence. Cranks gears have self-lock system. Hand cranks can be foldable and hidden while not in use. That prevents injuries from accidental and unwanted hits.

Mobilisation: If you need to move a patient with the bed, then It is a must to equip beds with castors. Especially for hospitals the mobility of beds is the basic need, yet for home care use we provide options for the beds to be without castors. Also more cost-effective option, beds without castors can be a good choice.

