

Thoracic Drainage

Ref: M4Z 200 0 0



For thoracic drainage procedures

- ▶ 3 L collection bottle made of polycarbonate material
- ▶ 25 cmH₂O vacuum power (at sea level)
- ▶ Intermittent vacuum
- ▶ Overflow safety device for pump
- ▶ Water manometer with cm scaled
- ▶ There is a safety jar that prevents fluid access to the vacuum pump from the water manometer
- ▶ Fan for cooling
- ▶ Phosphating process is applied before dyeing and powder paint (epoxy polyester) is housing
- ▶ Mobile
- ▶ Supply voltage: 230V, 50/60 Hz
- ▶ Energy consumption: 350 W
- ▶ Dimensions: 485 x 660 x 965 mm
- ▶ Weight: 28 kg
- ▶ Packing dimensions:
645 x 635 x 1030 mm

