

BIOBAG

Dialife

SWISS QUALITY DIALYSIS 



**SODIUM BICARBONATE BAG
FOR THE ON-LINE PRODUCTION
OF BICARBONATE CONCENTRATE**



INCREASED HYGIENE AND CLEANLINESS,
TO REDUCE RISK OF CONTAMINATION

Sodium bicarbonate bags are used for the on-line preparation of liquid bicarbonate concentrate via the hemodialysis machine. Treated water is channeled through the device under monitoring of the hemodialysis machine, to obtain a saturated bicarbonate solution.

Systems involving sodium bicarbonate bags represent a significant step forward for dialysis treatment, as they allow for a more hygienic supply of bicarbonate solution, hence reducing inflammations related to dialysis treatment while contributing to ultrapure dialysis.

Furthermore, the improved ergonomics of sodium bicarbonate bags offer additional convenience during usage and imply less storage space requirements.

Dialife offers high quality bicarbonate bags in the following sizes and formats:

PRODUCT CODE	CONTENT
BIOBAG DLB-650	650 gr. NaHCO3, Ph. Eur., USP
BIOBAG DLB-900	900 gr. NaHCO3, Ph. Eur., USP



LOWER SIZE AND WEIGHT,
FOR EASIER CARRIAGE, USE AND STORAGE

