

BIOCART

Dialife

SWISS QUALITY DIALYSIS 



**SODIUM BICARBONATE CARTRIDGE
FOR THE ON-LINE PRODUCTION OF
BICARBONATE CONCENTRATE**



 MANUFACTURED BY:
DIALIFE SA
Industrial Center Vedeggio 3
Via al Fiume 3
6807 Taverne
Switzerland


0476
UNI CEI EN ISO 13485:2016



www.dialifegroup.com

INCREASED HYGIENE AND CLEANLINESS, TO REDUCE RISK OF CONTAMINATION

Sodium bicarbonate cartridges are used for the on-line preparation of liquid bicarbonate concentrate via the hemodialysis machine. Treated water is channeled through the device under monitoring of the hemodialysis machine, to obtain a saturated bicarbonate solution.

Systems involving sodium bicarbonate cartridges represent a significant step forward for dialysis treatment, as they allow for a more hygienic supply of bicarbonate solution, hence reducing inflammations related to dialysis treatment while contributing to ultrapure dialysis.

Furthermore, the improved ergonomics of sodium bicarbonate cartridges offer additional convenience during usage and imply less storage space requirements.

Dialife offers high quality bicarbonate cartridges in the following sizes and formats:

PRODUCT CODE	CONTENT
BIOCART DLC-650	650 gr. NaHCO ₃ , Ph. Eur., USP
BIOCART DLC-720	720 gr. NaHCO ₃ , Ph. Eur., USP
BIOCART DLC-760	760 gr. NaHCO ₃ , Ph. Eur., USP
BIOCART DLC-920	920 gr. NaHCO ₃ , Ph. Eur., USP



LOWER SIZE AND WEIGHT, FOR EASIER CARRIAGE, USE AND STORAGE

