

DIADISINFECT

Dialife
SWISS QUALITY DIALYSIS 



DISINFECTANT FOR HEMODIALYSIS EQUIPMENT





DIADISINFECT 5L CANISTER



DIADISINFECT PLUS 10L CANISTER

Dialife disinfectants for hemodialysis machines are used to remove calcium/magnesium deposits and have fungicidal and bactericidal properties.

The product is pumped automatically by the hemodialysis device during the disinfection cycle.

DIADISINFECT contains citric acid, adjuvants and water.

DIADISINFECT PLUS contains citric acid, malic acid, lactic acid, adjuvants and water.

Both disinfectants provide decalcification at 37°C and disinfection starting at 60°C.

Available in 5 liters and 10 liters plastic canisters with the following compositions:

- DIADISINFECT-5 and DIADISINFECT-10: 50 % citric acid, adjuvants, water.
- DIADISINFECT plus-5 and DIADISINFECT plus-10: 21% citric acid, 2.5 % malic acid, 2.5 % lactic acid, adjuvants, water.

