


DIALUTION HD CONCENTRATE

Dialife
SWISS QUALITY DIALYSIS 



SODIUM BICARBONATE & ACID HEMODIALYSIS CONCENTRATE



 MANUFACTURED BY:
DIALIFE SA
Industrial Center Veduggio 3
Via al Fiume 3
6807 Taverne
Switzerland

CE
0476

UNI CEI EN ISO 13485:2016



www.dialifegroup.com



BICARBONATE CONCENTRATE



ACID CONCENTRATE

Bicarbonate and acid concentrates are mixed with purified water to generate dialysate solution, employed during hemodialysis treatment.

Our range of hemodialysis concentrates is European Pharmacopoeia compliant and undergoes strict production and quality control processes to ensure the highest safety and performance standards.

Dialife's DIALUTION solutions are offered as 5L or 10L canisters and with different formulation ranges.

With mixing ratios of 1/34 and 1/44, our acid concentrates are designed to be used with bicarbonate cartridges or liquid bicarbonate concentrate.

ACID CONCENTRATE FORMULATION RANGE (in mmol/L)

Sodium	138, 140
Potassium	0, 1, 1.5, 2, 3
Calcium	1, 1.25, 1.5, 1.75
Magnesium	0.5, 0.75, 1
Chloride	105.5 - 112.5
Acetic Acid	2, 3, 4, 5
Bicarbonate	32, 33, 34, 35, 36
Glucose	0, 5.55, 11.1

