

# Next Generation Digital Capillaroscopy!

Cost effective

Most popular

Plug & Play

All Inclusive

Capillaroscope  
Light

Capillaroscope  
Self-Install

Capillaroscope  
Pre-Installed

Capillaroscope  
Extensive

CAP-1

CAP-2

CAP-3

CAP-4




Item No.	Item Name
IN-007	Digital Capillaroscope, USB3.0
IN-010	20-50x Lens
IN-011	100x Lens
IN-012	200x Lens
IN-024	Immersion Oil
IN-025	Metris Software Light
IN-026	Metris Nailfold Capillaroscopy
IN-033	Single pedal Foot Switch
IN-034	Double pedal Foot Switch
IN-035	Capillaroscope XY-Stage
IN-036	Capillaroscope Stand
IN-037	Carrying Case
IN-038	Laptop or PC
IN-040	Healt care Cart


Non-invasive and safe method for morphologic examination and analysis of microcirculation abnormalities related to rheumatic disease, especially systemic sclerosis.

Metris Capillaroscopy software, offers the doctors intuitive yet powerful tools for early detection of digital ulcers.





Normal Pattern



Early Scleroderma Pattern



Active Scleroderma Pattern



Late Scleroderma Pattern



# Capillaroscope

Nail fold Video Capillaroscopy  
for Scleroderma Screening



V1, 2018. We reserve the right to modify data and specifications without prior notice.

Delivered by world's leading  
expert in Capillaroscopy  
Microscopes and Systems!