



Normal Pattern



Early Scleroderma Pattern



Active Scleroderma Pattern



Late Scleroderma Pattern

Images courtesy of Florian Meier and Inspectis.



Capillaroscope

Nailfold Video Capillaroscopy for
Scleroderma Screening




Powered by INSPECTIS® CAP Pro Software

Delivered by world's leading expert in
Capillaroscopy Microscopes and Systems!

V1, 2020. We reserve the right to modify data and specifications without prior notice.

New Generation Digital Capillaroscopy!



Cost effective	Most popular	Plug & Play	All Inclusive	
Capillaroscope Light	Capillaroscope Self-Install	Capillaroscope Pre-Installed	Capillaroscope Extensive	
CAP-1	CAP-2	CAP-3	CAP-4	
Item No.	Item Name			
IN-001	Digital Capillaroscope, USB3.0	●	●	●
IN-012	200x Lens	●	●	●
IN-024	Immersion Oil	●	●	●
IN-025	INSPECTIS® CAP Basics software	●		
IN-026	INSPECTIS® CAP Pro software		●	●
IN-033	Single-pedal Foot Switch	●		
IN-034	Triple-pedal Foot Switch		●	●
IN-036	Capillaroscope bench-top Holder	●	●	●
IN-037	Carrying Case		●	●
IN-038	Laptop or PC			●
IN-039	Laptop soft case			●
IN-040	Health care Cart			●
IN-010	15-45x Varifocal Lens with Polariser	Optional		Optional
IN-011	100x Lens	Optional		Optional
IN-035	Capillaroscope Stand with XY-stage	Optional		Optional

Non-invasive and safe method for morphologic examination and analysis of microcirculation abnormalities related to rheumatic disease, especially systemic sclerosis.

INSPECTIS® CAP Pro software offers the doctors intuitive yet powerful tools for efficient Nailfold capillaroscopy.

